

MED30 recombinant protein

Recombinant Human MED30

CATALOG #
AAA250787

HOST
E Coli


PRODUCT SYNONYMS

MED30; Recombinant Human MED30; mediator complex subunit 30; THRAP6; TRAP25; MED30 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli
Purity/Purification	Greater than 85% as determined by SDS-PAGE
Form/Format	Lyophilized from a 0.2 um filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH7.4
Sequence Positions	2-178
Species	Human
Tag	N-terminal His-IF2DI Tag
Reconstitution	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to a concentration of 0.1-1 mg/ml by gently pipetting 2-3 times, do not vortex
Preparation and Storage	The lyophilized protein is stable at -20°C for up to 1 year. The protein solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycles.

IMAGES



SDS-PAGE

NCBI AND UNIPROT INFORMATION

NCBI GI #	545479145
NCBI GeneID	90390
NCBI Accession #	NP_001269915.1

NCBI GenBank Nucleotide #	NM_001282986.1
UniProt Accession #	Q96HR3
Molecular Weight	40.8 kDa under reducing conditions
NCBI Official Full Name	mediator of RNA polymerase II transcription subunit 30 isoform 2
NCBI Official Synonym Full Names	mediator complex subunit 30
NCBI Official Symbol	MED30
NCBI Official Synonym Symbols	MED30S; THRAP6; TRAP25
NCBI Protein Information	mediator of RNA polymerase II transcription subunit 30; TRAP/Mediator complex component TRAP25; thyroid hormone receptor associated protein 6; thyroid hormone receptor-associated protein 6; putative mediator of RNA polymerase II transcription subunit 30;
UniProt Protein Name	Mediator of RNA polymerase II transcription subunit 30
UniProt Gene Name	MED30
UniProt Synonym Gene Names	THRAP6; TRAP25; Trap25
UniProt Entry Name	MED30_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/med30-recombinant-protein-e-coli/250787>