

Human COVID-19 Nucleocapsid (NP) Spike RBD Alpha B.1.1.7 Variant (UK Variant) IgM Coronavirus ELISA Kit | COVID-19 elisa kit

Human anti-SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant (UK Variant) IgM ELISA Kit

CATALOG #
AAA250912

REACTIVITY
Human


PRODUCT SYNONYMS

COVID-19 Nucleocapsid (NP) Spike RBD Alpha B.1.1.7 Variant (UK Variant) IgM Coronavirus; Human anti-SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant (UK Variant) IgM ELISA Kit; 2019 Novel Coronavirus; Coronavirus; CoV; COVID-19 virus; HCoV-2; Human Coronavirus 2019; SARS2; SARS-CoV-2; Severe acute respiratory syndrome coronavirus 2; anti-SARS-CoV-2 UK Variant IgM; Spike RBD Alpha B.1.1.7 Variant (UK Variant) IgM; COVID-19 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Anti SARS-CoV-2 UK Variant IgM. No significant cross-reactivity or interference between Anti SARS-CoV-2 UK Variant IgM and analogues was observed.
Samples	Serum, plasma and other biological fluids
Detection Range	0.78-50ng/ml
Sensitivity	0.469ng/ml
Intra-assay Precision	CV<8%
Inter-assay Precision	CV<10%
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for COVID-19 elisa kit	Principle of the Assay: This kit was based on indirect enzyme-linked immune-sorbent assay technology. Recombinant nucleocapsid protein was precoated onto 96-well plates. And the HRP conjugated antibody was used as detection antibodies. The standards, test samples and HRP conjugated detection antibody were added to the wells subsequently, and washed with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the target amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of target can be calculated.
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/covid-19-nucleocapsid-np-spike-rbd-alpha-b-1-1-7-variant-uk-variant-igm-coronavirus/250912>