

# Human TCC-CH50 (Terminal Complement Complex-50% complement hemolysis) ELISA Kit | TCC-CH50 elisa kit

Human TCC-CH50 (Terminal Complement Complex-50% complement hemolysis) ELISA Kit

CATALOG #  
**AAA250982**

REACTIVITY  
**Human**


## PRODUCT SYNONYMS

TCC-CH50 (Terminal Complement Complex-50% complement hemolysis); Human TCC-CH50 (Terminal Complement Complex-50% complement hemolysis) ELISA Kit; TCC-CH50; Terminal Complement Complex-50% complement hemolysis; TCC-CH50 elisa kit

## PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                              |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Human                                                                                                                                                        |
| Specificity             | Specifically recognize CH50, no obvious cross reaction with other analogues                                                                                  |
| Assay Type              | Sandwich                                                                                                                                                     |
| Samples                 | Serum, plasma, cell culture supernatant and other biological samples                                                                                         |
| Detection Range         | 4.69-300U/ml                                                                                                                                                 |
| Sensitivity             | 2.813U/ml                                                                                                                                                    |
| Intra-assay Precision   | Intra-assay Precision: samples with low, medium and high concentration are tested 20 times on same plate.                                                    |
| Inter-assay Precision   | Inter-assay Precision: samples with low, medium and high concentration are tested 20 times on three different plates.                                        |
| Preparation and Storage | Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C |

## IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|----------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for TCC-CH50 elisa kit | Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti CH50 antibody was pre-coated onto the 96-well plate. The biotin conjugated anti CH50 antibody was used as the detection antibody. The standards and pilot samples were added to the wells subsequently. After incubation, unbound conjugates were removed by wash buffer. Then, biotinylated detection antibody was added to bind with CH50 conjugated on coated antibody. After washing off unbound conjugates, HRP-Streptavidin was added. After a third washing, TMB substrates were added to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that turned yellow after adding a stop solution. Read the O.D. absorbance at 450nm in a microplate reader. The concentration of CH50 in the sample was calculated by drawing a standard curve. The concentration of the target substance is proportional to the OD450 value. |
|----------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tcc-ch50-terminal-complement-complex-50-complement-hemolysis/250982>