

Mouse DHEA-S (Dehydroepiandrosterone sulfate) ELISA Kit | DHEA-S elisa kit

Mouse DHEA-S (Dehydroepiandrosterone sulfate) ELISA Kit

CATALOG #
AAA250993

REACTIVITY
Mouse


PRODUCT SYNONYMS

DHEA-S (Dehydroepiandrosterone sulfate); Mouse DHEA-S (Dehydroepiandrosterone sulfate) ELISA Kit; DHEA-S; DHEAS; DHEA-ST; hSTa; ST2A3; SULT2A1; alcohol; hydroxysteroid sulfotransferase; Dehydroepiandrosterone sulfotransferase; DHEA-STbile salt sulfotransferase; EC 2.8.2.14; HST; Hydroxysteroid Sulfotransferase; ST2; ST2A1; STDbile-salt sulfotransferase 2A1; Sulfotransferase 2A1; sulfotransferase family; cytosolic; 2A; dehydroepiandrosterone; DHEA-preferring; member 1; DHEA-S elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	Specifically recognize DHEA-S, no obvious cross reaction with other analogues
Assay Type	Competitive
Samples	Serum, plasma, cell culture supernatant and other biological samples
Detection Range	1.563-100ng/ml
Sensitivity	0.938ng/ml
Intra-assay Precision	Intra-assay Precision: samples with low, medium and high concentration are tested 20 times on same plate.
Inter-assay Precision	Inter-assay Precision: samples with low, medium and high concentration are tested 20 times on three different plates.
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DHEA-S elisa kit	Principle of the Assay: This kit was based on Competitive-ELISA detection method. The microtiter plate provided in this kit has been pre-coated with DHEA-S. During the reaction, DHEA-S in the sample or standard competes with a fixed amount of DHEA-S on the solid phase supporter for sites on the Biotinylated Detection Antibody specific to DHEA-S. Excess conjugate and unbound sample or standard are washed from the plate, and HRP-Streptavidin (SABC) is added to each microplate well and incubated. Then TMB substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of a acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm. The concentration of DHEA-S in the samples is then determined by comparing the OD of the samples to the standard curve. The concentration of the target substance was inversely proportional to the OD450 value.
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/dhea-s-dehydroepiandrosterone-sulfate/250993>