

Human TIRAP (Toll/interleukin-1 receptor domain-containing adapter protein) ELISA Kit | TIRAP elisa kit

Human TIRAP (Toll/interleukin-1 receptor domain-containing adapter protein) QuickTest ELISA Kit

CATALOG # AAA252563	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

TIRAP (Toll/interleukin-1 receptor domain-containing adapter protein); Human TIRAP (Toll/interleukin-1 receptor domain-containing adapter protein) QuickTest ELISA Kit; Toll/interleukin-1 receptor domain-containing adapter protein; TIR domain-containing adapter protein; Adaptor protein Wyatt; MyD88 adapter-like protein; MyD88-2; TIRAP; MAL; TIRAP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Detection Range	31.25-2000pg/ml
Sensitivity	18.75pg/ml
Time	120 minutes
Preparation and Storage	Store at 2-8 degree C for 12 months

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	89111122
NCBI GeneID	114609
NCBI Accession #	NP_001034750.1
NCBI GenBank Nucleotide #	NM_001039661.1
UniProt Accession #	P58753 ; Q56UH9 ; Q56UI0 ; Q8N5E5 ; B3KW65
Molecular Weight	23,883 Da
NCBI Official Full Name	toll/interleukin-1 receptor domain-containing adapter protein isoform a
NCBI Official Synonym Full Names	toll-interleukin 1 receptor (TIR) domain containing adaptor protein
NCBI Official Symbol	TIRAP
NCBI Official Synonym Symbols	Mal; wyatt; BACTS1; MyD88-2

NCBI Protein Information	toll/interleukin-1 receptor domain-containing adapter protein; adapter protein wyatt; adaptor protein Wyatt; MyD88 adapter-like protein; Toll-like receptor adaptor protein
UniProt Protein Name	Toll/interleukin-1 receptor domain-containing adapter protein
UniProt Gene Name	TIRAP
UniProt Synonym Gene Names	MAL; TIR domain-containing adapter protein; MyD88-2
UniProt Entry Name	TIRAP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tirap-toll-interleukin-1-receptor-domain-containing-adapter-protein/252563>