

Human PIP4K2C (Phosphatidylinositol-5-phosphate 4-kinase type-2 gamma) ELISA Kit | PIP4K2C elisa kit

Human PIP4K2C (Phosphatidylinositol-5-phosphate 4-kinase type-2 gamma) QuickTest ELISA Kit

CATALOG # AAA252566	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

PIP4K2C (Phosphatidylinositol-5-phosphate 4-kinase type-2 gamma); Human PIP4K2C (Phosphatidylinositol-5-phosphate 4-kinase type-2 gamma) QuickTest ELISA Kit; Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma; Phosphatidylinositol 5-phosphate 4-kinase type II gamma; PI(5)P 4-kinase type II gamma; PIP4KII-gamma; PIP4K2C; PIP5K2C; PIP4K2C elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Detection Range	78.125-5000pg/ml
Sensitivity	46.875pg/ml
Time	120 minutes
Preparation and Storage	Store at 2-8 degree C for 12 months

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	226371739
NCBI GeneID	79837
NCBI Accession #	NP_001139730
NCBI GenBank Nucleotide #	NM_001146258
UniProt Accession #	Q8TBX8 ; Q9H6N2 ; B2RDL3 ; B4DM11 ; B4DY44
Molecular Weight	45,053 Da
NCBI Official Full Name	phosphatidylinositol 5-phosphate 4-kinase type-2 gamma isoform a
NCBI Official Synonym Full Names	phosphatidylinositol-5-phosphate 4-kinase, type II, gamma
NCBI Official Symbol	PIP4K2C
NCBI Official Synonym Symbols	PIP5K2C

NCBI Protein Information	phosphatidylinositol 5-phosphate 4-kinase type-2 gamma; PI(5)P 4-kinase type II gamma; PIP4KII-gamma; phosphatidylinositol 4-phosphate 5-kinase; phosphatidylinositol 5-phosphate 4-kinase type II gamma; phosphatidylinositol-4-phosphate 5-kinase, type II, g
UniProt Protein Name	Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma
UniProt Gene Name	PIP4K2C
UniProt Synonym Gene Names	PIP5K2C; PI(5)P 4-kinase type II gamma; PIP4KII-gamma
UniProt Entry Name	PI42C_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/pip4k2c-phosphatidylinositol-5-phosphate-4-kinase-type-2-gamma/252566>