

Human FCER1A (High affinity immunoglobulin epsilon receptor subunit alpha) ELISA Kit | FCER1A elisa kit

Human FCER1A (High affinity immunoglobulin epsilon receptor subunit alpha) QuickTest ELISA Kit

CATALOG # AAA252568	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

FCER1A (High affinity immunoglobulin epsilon receptor subunit alpha); Human FCER1A (High affinity immunoglobulin epsilon receptor subunit alpha) QuickTest ELISA Kit; High affinity immunoglobulin epsilon receptor subunit alpha; Fc-epsilon RI-alpha; FcERI; IgE Fc receptor subunit alpha; FCER1A; FCE1A; FCER1A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Detection Range	31.25-2000pg/ml
Sensitivity	18.75pg/ml
Time	120 minutes
Preparation and Storage	Store at 2-8 degree C for 12 months

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	119865
NCBI GeneID	2205
UniProt Accession #	P12319
Molecular Weight	30kDa
NCBI Official Full Name	High affinity immunoglobulin epsilon receptor subunit alpha
NCBI Official Synonym Full Names	Fc fragment of IgE, high affinity I, receptor for; alpha polypeptide

NCBI Official Symbol	FCER1A
NCBI Official Synonym Symbols	FCE1A; FcERI
NCBI Protein Information	high affinity immunoglobulin epsilon receptor subunit alpha; Fc-epsilon RI-alpha; Fc epsilon RI alpha-chain; IgE Fc receptor subunit alpha; Fc IgE receptor, alpha polypeptide; high affinity immunoglobulin epsilon receptor alpha-subunit; immunoglobulin E r
UniProt Protein Name	High affinity immunoglobulin epsilon receptor subunit alpha
UniProt Gene Name	FCER1A
UniProt Synonym Gene Names	FCE1A; FcERI
UniProt Entry Name	FCERA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/fcer1a-high-affinity-immunoglobulin-epsilon-receptor-subunit-alpha/252568>