

Human PARP2 (Poly [ADP-ribose] polymerase 2) ELISA Kit | PARP2 elisa kit

Human PARP2 (Poly [ADP-ribose] polymerase 2) QuickTest ELISA Kit

CATALOG #
AAA252574

REACTIVITY
Human


PRODUCT SYNONYMS

PARP2 (Poly [ADP-ribose] polymerase 2); Human PARP2 (Poly [ADP-ribose] polymerase 2) QuickTest ELISA Kit; Poly [ADP-ribose] polymerase 2; PARP-2; hPARP-2; ADP-ribosyltransferase diphtheria toxin-like 2; ARTD2; DNA ADP-ribosyltransferase PARP2; NAD(+) ADP-ribosyltransferase 2; ADPRT-2; Poly[ADP-ribose] synthase 2; pADPRT-2; Protein poly-ADP-ribosyltransferase PARP2; PARP2; ADPRT2; ADPRTL2; PARP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Detection Range	0.313-20ng/ml
Sensitivity	0.188ng/ml
Time	120 minutes
Preparation and Storage	Store at 2-8 degree C for 12 months

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	110825963
NCBI GeneID	10038
NCBI Accession #	NP_001036083.1
NCBI GenBank Nucleotide #	NM_001042618.1
UniProt Accession #	Q9UGN5 ; Q8TEU4 ; Q9NUV2 ; Q9UMR4 ; Q9Y6C8
Molecular Weight	66,206 Da
NCBI Official Full Name	poly

NCBI Official Synonym Full Names	poly (ADP-ribose) polymerase 2
NCBI Official Symbol	PARP2
NCBI Official Synonym Symbols	ARTD2; ADPRT2; PARP-2; ADPRTL2; ADPRTL3; pADPRT-2
NCBI Protein Information	poly [ADP-ribose] polymerase 2; poly [ADP-ribose] polymerase 2; ADPRT-2; hPARP-2; poly(ADP-ribose) synthetase; poly[ADP-ribose] synthase 2; poly[ADP-ribose] synthetase 2; NAD(+) ADP-ribosyltransferase 2; poly (ADP-ribosyl) transferase-like 2; poly (ADP-ri
UniProt Protein Name	Poly [ADP-ribose] polymerase 2
UniProt Gene Name	PARP2
UniProt Synonym Gene Names	ADPRT2; ADPRTL2; PARP-2; hPARP-2; ARTD2; ADPRT-2; pADPRT-2
UniProt Entry Name	PARP2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/parp2-poly-adp-ribose-polymerase-2/252574>