

Human HLA-DQB2 (HLA class II histocompatibility antigen, DQ beta 2 chain) ELISA Kit | HLA-DQB2 elisa kit

Human HLA-DQB2 (HLA class II histocompatibility antigen, DQ beta 2 chain) QuickTest ELISA Kit

CATALOG #
AAA252626

REACTIVITY
Human


PRODUCT SYNONYMS

HLA-DQB2 (HLA class II histocompatibility antigen; DQ beta 2 chain); Human HLA-DQB2 (HLA class II histocompatibility antigen; DQ beta 2 chain) QuickTest ELISA Kit; HLA class II histocompatibility antigen; DQ beta 2 chain; DX beta chain; MHC class II antigen DQB2; HLA-DQB2; HLA-DXB; HLA-DQB2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Detection Range	0.313-20ng/ml
Sensitivity	0.188ng/ml
Time	120 minutes
Preparation and Storage	Store at 2-8 degree C for 12 months

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	311893321
NCBI GeneID	3120
NCBI Accession #	NP_001185787.1
NCBI GenBank Nucleotide #	NM_001198858.1
UniProt Accession #	P05538; Q29826; Q29870; Q29871; Q29872; Q29873; Q5SR06; A6NIA5
Molecular Weight	26,031 Da

NCBI Official Full Name	HLA class II histocompatibility antigen, DQ beta 2 chain isoform 2
NCBI Official Synonym Full Names	major histocompatibility complex, class II, DQ beta 2
NCBI Official Symbol	HLA-DQB2
NCBI Official Synonym Symbols	HLA-DXB; HLA-DQB1
NCBI Protein Information	HLA class II histocompatibility antigen, DQ beta 2 chain; MHC class II antigen DQB2; major histocompatibility complex, class II, DQ beta 1; HLA class II histocompatibility antigen, DX beta chain; DV19.1 (major histocompatibility complex, class II, DQ beta
UniProt Protein Name	HLA class II histocompatibility antigen, DQ beta 2 chain
UniProt Gene Name	HLA-DQB2
UniProt Synonym Gene Names	HLA-DXB
UniProt Entry Name	DQB2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/hla-dqb2-hla-class-ii-histocompatibility-antigen-dq-beta-2-chain/252626>