

## Rat INSR-beta (Insulin Receptor Beta Subunit) ELISA Kit | INSR-beta elisa kit

Rat INSR-beta (Insulin Receptor Beta Subunit) QuickTest ELISA Kit

CATALOG #  
**AAA252863**

REACTIVITY  
**Rat**


### PRODUCT SYNONYMS

INSR-beta (Insulin Receptor Beta Subunit); Rat INSR-beta (Insulin Receptor Beta Subunit) QuickTest ELISA Kit; INSR-beta; Insulin Receptor Beta Subunit; CD220; Insr; INSR-beta elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                       |
|-------------------------|-------------------------------------------------------------------------------------------------------|
| Reactivity              | Rat                                                                                                   |
| Assay Type              | Sandwich ELISA, Double Antibody                                                                       |
| Samples                 | Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples |
| Detection Range         | 0.938-60ng/ml                                                                                         |
| Sensitivity             | 0.563ng/ml                                                                                            |
| Time                    | 120 minutes                                                                                           |
| Preparation and Storage | Store at 2-8 degree C for 12 months                                                                   |

### IMAGES



Standard Curve (Sample)

### NCBI AND UNIPROT INFORMATION

|                            |                        |
|----------------------------|------------------------|
| UniProt Accession #        | <a href="#">P15127</a> |
| UniProt Protein Name       | Insulin receptor       |
| UniProt Gene Name          | <b>Insr</b>            |
| UniProt Synonym Gene Names | IR                     |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.