

Human PSME2 (Proteasome activator complex subunit 2) ELISA Kit | PSME2 elisa kit

Human PSME2 (Proteasome activator complex subunit 2) ELISA Kit

CATALOG #
AAA253689

REACTIVITY
Human


PRODUCT SYNONYMS

PSME2 (Proteasome activator complex subunit 2); Human PSME2 (Proteasome activator complex subunit 2) ELISA Kit; Proteasome activator complex subunit 2; 11S regulator complex subunit beta; REG-beta; Activator of multicatalytic protease subunit 2; Proteasome activator 28 subunit beta; PA28b; PA28beta; PSME2; PSME2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	31.25-2000pg/ml
Sensitivity	18.75pg/ml
Time	4 hours
Preparation and Storage	Store at 2-8 degree C for 12 months

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	30410792
NCBI GeneID	5721
NCBI Accession #	NP_002809.2
NCBI GenBank Nucleotide #	NM_002818.2
UniProt Accession #	Q9UL46; Q15129
Molecular Weight	27,402 Da
NCBI Official Full Name	proteasome activator complex subunit 2

NCBI Official Synonym Full Names	proteasome (prosome, macropain) activator subunit 2 (PA28 beta)
NCBI Official Symbol	PSME2
NCBI Official Synonym Symbols	PA28B; REGbeta; PA28beta
NCBI Protein Information	proteasome activator complex subunit 2; 11S regulator complex beta subunit; 11S regulator complex subunit beta; MCP activator, 31-kD subunit; REG-beta; activator of multicatalytic protease subunit 2; cell migration-inducing protein 22; proteasome activato
UniProt Protein Name	Proteasome activator complex subunit 2
UniProt Gene Name	PSME2
UniProt Synonym Gene Names	REG-beta; PA28b; PA28beta
UniProt Entry Name	PSME2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/psme2-proteasome-activator-complex-subunit-2/253689>