

Rabbit anti-Human AMIGO1 Polyclonal Antibody | anti-AMIGO1 antibody

Anti-AMIGO1 Antibody, Rabbit Polyclonal


CATALOG # AAA258975	HOST Rabbit	REACTIVITY Human	APPLICATIONS Immunohistochemistry, Western Blot
-------------------------------	-----------------------	----------------------------	---

PRODUCT SYNONYMS


AMIGO1; Polyclonal Antibody; Anti-AMIGO1 Antibody; Rabbit Polyclonal; adhesion molecule with Ig-like domain 1; Anti-ALI2 Antibody; Anti-AMIGO Antibody; Anti-AMIGO-1 Antibody; anti-AMIGO1 antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	Human
Clonality	Polyclonal
Isotype	Rabbit IgG
Specificity	Human AMIGO1
Purity/Purification	Protein A & Antigen Affinity
Form/Format	Liquid; PBS, pH7.0 with 0.03% Proclin300
Applicable Applications for anti-AMIGO1 antibody	IHC (Immunohistochemistry), WB (Western Blot)
Immunogen	E coli-derived Human AMIGO1 fragment
Preparation	Produced in rabbits immunized with E coli-derived Human AMIGO1 fragment, and purified by antigen affinity chromatography.
Preparation and Storage	This antibody can be stored at 2-8 degree C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20 degree C to -80 degree C. Avoid repeated freeze-thaw cycles.



IHC (Immunohistochemistry)
Immunochemical staining of human AMIGO1 in human breast with rabbit polyclonal antibody at 1:500 dilution, formalin-fixed paraffin embedded sections.



WB (Western Blot)
Anti-AMIGO1 rabbit polyclonal antibody at 1:500 dilutionLane A: Hela Whole Cell LysateLysates/proteins at 30 ug per lane.SecondaryGoat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.Developed using the ECL technique.Performed under reducing conditions.Predicted band size:55 kDaObserved band size:55 kDa(We are unsure as to the identity of these extra bands.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/polyclonal/amigo1-antibody-immunohistochemistry-western-blot-elisa/258975>