

Rabbit anti-Human PFDN5 Polyclonal Antibody | anti-PFDN5 antibody

Anti-PFDN5 Antibody, Rabbit Polyclonal


CATALOG # AAA259071	HOST Rabbit	REACTIVITY Human	APPLICATIONS Immunoprecipitation, Western Blot
-------------------------------	-----------------------	----------------------------	--

PRODUCT SYNONYMS


PFDN5; Polyclonal Antibody; Anti-PFDN5 Antibody; Rabbit Polyclonal; prefoldin subunit 5; Anti-MM-1 Antibody; Anti-MM1 Antibody; Anti-PFD5 Antibody; anti-PFDN5 antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	Human
Clonality	Polyclonal
Isotype	Rabbit IgG
Specificity	Human PFDN5
Purity/Purification	Protein A & Antigen Affinity
Form/Format	Liquid; PBS, pH7.0 with 0.03% Proclin300
Applicable Applications for anti-PFDN5 antibody	IP (Immunoprecipitation), WB (Western Blot)
Immunogen	E coli-derived Human PFDN5 fragment
Preparation	Produced in rabbits immunized with E coli-derived Human PFDN5 fragment, and purified by antigen affinity chromatography.
Preparation and Storage	This antibody can be stored at 2-8 degree C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20 degree C to -80 degree C. Avoid repeated freeze-thaw cycles.



IP (Immunoprecipitation)
PFDN5 was immunoprecipitated using: Lane A: 0.5 mg HepG2 Whole Cell Lysate 4 uL anti-PFDN5 rabbit polyclonal antibody and 60 ug of Immunomagnetic beads Protein A/G. Primary antibody: Anti-PFDN5 rabbit polyclonal antibody, at 1:100 dilution. Secondary antibody: Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. Developed using the ECL technique. Performed under reducing conditions. Predicted band size: 17 kDa. Observed band size: 15 kDa.



WB (Western Blot)
Anti-PFDN5 rabbit polyclonal antibody at 1:500 dilution. Lane A: HeLa Whole Cell Lysate. Lane B: HepG2 Whole Cell Lysate. Lane C: 293T Whole Cell Lysate. Lane D: Jurkat Whole Cell Lysate. Lysates/proteins at 30 ug per lane. Secondary: Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. Developed using the ECL technique. Performed under reducing conditions. Predicted band size: 17 kDa. Observed band size: 15 kDa. (We are unsure as to the identity of these extra bands.)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.