

General GAPDH Assay Kit | GAPDH assay kit

GAPDH Activity Microplate Assay Kit

CATALOG # AAA266483	REACTIVITY General	APPLICATIONS Functional Assay
-------------------------------	------------------------------	---


PRODUCT SYNONYMS

GAPDH; GAPDH Activity Microplate Assay Kit; EC 1.2.1.12; glyceraldehyde-3-phosphate dehydrogenase; NAD-dependent glyceraldehyde-3-phosphate dehydrogenase; triosephosphate dehydrogenase; GAPDH assay kit

PRODUCT SPECIFICATIONS


Reactivity	General
Samples	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Detection Method	Colorimetric
Detection Range	4umol/L-400umol/L
CAS Number	9001-50-7
Platform	Microplate reader 492nm
Applicable Applications for GAPDH assay kit	FA (Functional Assay)
Preparation and Storage	Shipped and store at 4 degree C for 6 months.

IMAGES



Application Data

Positive Control reaction in 96-well plate assay with decreasing the concentration



Standard Curve (Sample)

Representative Standard Curve

ADDITIONAL INFORMATION

Related Product Information for GAPDH assay kit	Detection and Quantification of GAPDH Activity.
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.