

General Trehalose Assay Kit

Trehalose Fluorometric Microplate Assay Kit

CATALOG # AAA266504	REACTIVITY General	APPLICATIONS Functional Assay
-------------------------------	------------------------------	---


PRODUCT SYNONYMS

Trehalose; Trehalose Fluorometric Microplate Assay Kit; mycose; tremalose; Trehalose assay kit

PRODUCT SPECIFICATIONS

Reactivity	General
Samples	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Detection Method	Fluorimetric
Detection Range	2umol/L-200umol/L
CAS Number	99-20-7
Platform	Fluorescence microplate reader Ex/Em = 535/587
Applicable Applications for Trehalose assay kit	FA (Functional Assay)
Preparation and Storage	Shipped and store at 4 degree C for 6 months.

IMAGES



Standard Curve (Sample)

Representative Standard Curve

ADDITIONAL INFORMATION

Related Product Information for Trehalose assay kit	Detection and Quantification of Trehalose Content.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/general-assay-kits/trehalose/266504>