

General Trehalase Activity Assay Kit

Trehalase Activity Fluorometric Microplate Assay Kit

CATALOG # AAA266505	REACTIVITY General	APPLICATIONS Functional Assay
-------------------------------	------------------------------	---


PRODUCT SYNONYMS

Trehalase Activity; Trehalase Activity Fluorometric Microplate Assay Kit; EC 3.2.1.28; Trehalase; Alpha; alpha-trehalase; TREH; Trehalase Activity assay kit

PRODUCT SPECIFICATIONS


Reactivity	General
Samples	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Detection Method	Fluorimetric
Detection Range	1umol/L-100umol/L
CAS Number	9025-52-9
Platform	Fluorescence microplate reader Ex/Em = 535/587
Applicable Applications for Trehalase Activity assay kit	FA (Functional Assay)
Preparation and Storage	Shipped and store at 4 degree C for 6 months.

IMAGES



Application Data

Positive Control reaction in 96-well plate assay with decreasing the concentration



Standard Curve (Sample)

Representative Standard Curve

ADDITIONAL INFORMATION

Related Product Information for Trehalase Activity assay kit	Detection and Quantification of Trehalase Activity.
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/general-assay-kits/trehalase-activity/266505>