

General Glucose Oxidase Activity Assay Kit

Glucose Oxidase Activity Fluorometric Microplate Assay Kit

CATALOG # AAA266508	REACTIVITY General	APPLICATIONS Functional Assay
-------------------------------	------------------------------	---


PRODUCT SYNONYMS

Glucose Oxidase Activity; Glucose Oxidase Activity Fluorometric Microplate Assay Kit; EC 1.1.3.4; Glucose oxidase; Beta-D-glucose:oxygen 1-oxido-reductase; D-glucose-1-oxidase; Glucose aerodehydrogenase; Glucose oxyhydrase; GOD; Glucose Oxidase Activity assay kit

PRODUCT SPECIFICATIONS

Reactivity	General
Samples	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Detection Method	Fluorimetric
Detection Range	1umol/L-100umol/L
CAS Number	9001-37-0
Platform	Fluorescence microplate reader Ex/Em = 535/587
Applicable Applications for Glucose Oxidase Activity assay kit	FA (Functional Assay)
Preparation and Storage	Shipped and store at 4 degree C for 6 months.

IMAGES



Standard Curve (Sample)
Representative Standard Curve

ADDITIONAL INFORMATION

Related Product Information for Glucose Oxidase Activity assay kit	Detection and Quantification of Glucose Oxidase Activity.
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.