

Pig HSPA1A (Heat Shock 70kDa Protein 1A) ELISA Kit | HSPA1A elisa kit

Pig HSPA1A (Heat Shock 70kDa Protein 1A) ELISA Kit

CATALOG #
AAA270165

REACTIVITY
Pig


PRODUCT SYNONYMS

HSPA1A (Heat Shock 70kDa Protein 1A); Pig HSPA1A (Heat Shock 70kDa Protein 1A) ELISA Kit; HSP70-1; HSP-A1A; HSPA1; HSP70-2; HSX70; HSPA1A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Pig
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Assay Type	Sandwich
Detection Range	1.57-100 ng/mL
Sensitivity	0.61 ng/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HSPA1A elisa kit

Principle of the assay: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Heat Shock 70kDa Protein 1A(HSPA1A). Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Heat Shock 70kDa Protein 1A(HSPA1A). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Heat Shock 70kDa Protein 1A(HSPA1A), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Heat Shock 70kDa Protein 1A(HSPA1A) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for HSPA1A elisa kit

Signal transduction; Tumor immunity; Infection immunity; Developmental science

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/pig-elisa-kits/hspa1a-heat-shock-70kda-protein-1a/270165>