

General gABA (Gamma-Aminobutyric Acid) ELISA Kit | gABA elisa kit

General gABA (Gamma-Aminobutyric Acid) ELISA Kit

CATALOG #
AAA270177

REACTIVITY
General


PRODUCT SYNONYMS

gABA (Gamma-Aminobutyric Acid); General gABA (Gamma-Aminobutyric Acid) ELISA Kit; 4-Aminobutyric Acid; gABA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	General
Samples	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Assay Type	Competitive
Detection Range	31.25-2000 pg/mL
Sensitivity	9.27 pg/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for gABA elisa kit

Principle of the assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Gamma-Aminobutyric Acid(gABA) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Gamma-Aminobutyric Acid(gABA). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Gamma-Aminobutyric Acid(gABA) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for gABA elisa kit

[Signal transduction](#); [Cardiovascular biology](#); [Neuro science](#); [Bone metabolism](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

