

Mouse ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit | ADAM8 elisa kit

Mouse ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit

CATALOG #
AAA270316

REACTIVITY
Mouse


PRODUCT SYNONYMS

ADAM8 (A Disintegrin And Metalloprotease 8); Mouse ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit; CD156-A; MS2; CD156a; Disintegrin and metalloproteinase domain-containing protein 8; ADAM8 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, tissue homogenates and other biological fluids
Assay Type	Sandwich
Detection Range	7.82-500 pg/mL
Sensitivity	3.3 pg/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ADAM8 elisa kit

Principle of the assay: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to A Disintegrin And Metalloprotease 8(ADAM8). Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to A Disintegrin And Metalloprotease 8(ADAM8). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain A Disintegrin And Metalloprotease 8(ADAM8), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of A Disintegrin And Metalloprotease 8(ADAM8) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for ADAM8 elisa kit

Enzyme & Kinase; Tumor immunity; Immune molecule; Neuro science

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/adam8-a-disintegrin-and-metalloprotease-8/270316>