

Human PCII (Procollagen II) ELISA Kit | PCII elisa kit

Human PCII (Procollagen II) ELISA Kit

CATALOG #
AAA271158

REACTIVITY
Human


PRODUCT SYNONYMS

PCII (Procollagen II); Human PCII (Procollagen II) ELISA Kit; PCOL-II; PC-II; PC2; Type II Procollagen; PCII elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma and other biological fluids.
Assay Type	Sandwich
Detection Range	0.32-20 ng/mL
Sensitivity	0.107 ng/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PCII elisa kit

Principle of the assay: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Procollagen II(PCII). Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Procollagen II(PCII). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Procollagen II(PCII), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Procollagen II(PCII) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for PCII elisa kit

Metabolic pathway; Infection immunity; Developmental science; Hepatology; Nutrition metabolism; Dermatology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/pcii-procollagen-iii/271158>