

Porcine Creatine Kinase MM Isoenzyme (CK MM) ELISA Kit | CK MM elisa kit

Porcine Creatine Kinase MM Isoenzyme (CK MM) ELISA Kit

CATALOG #
AAA27439

REACTIVITY
Porcine


PRODUCT SYNONYMS

Creatine Kinase MM Isoenzyme (CK MM); Porcine Creatine Kinase MM Isoenzyme (CK MM) ELISA Kit; CK MM elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Specificity	This assay has high sensitivity and excellent specificity for detection of CK-MM. No significant cross-reactivity or interference between CK-MM and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between CK-MM and all the analogues, therefore, cross reaction may still exist in some cases.
Assay Type	Quantitative Competitive
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Sensitivity	1.0 ng/mL
Reactivity Note	No significant cross-reactivity or interference between this antigen and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection, therefore, cross reaction may still exist in some cases.
Preparation and Storage	Store at 2-8 degree C. The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CK MM elisa kit	Principle of the Assay: CK-MM ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-CK-MM antibody and an CK-MM-HRP conjugate. The assay sample and buffer are incubated together with CK-MM-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the CK-MM concentration since CK-MM from samples and CK-MM-HRP conjugate compete for the anti-CK-MM antibody binding site. Since the number of sites is limited, as more sites are occupied by CK-MM from the sample, fewer sites are left to bind CK-MM-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The CK-MM concentration in each sample is interpolated from this standard curve. Intended Uses: This CK-MM ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Porcine CK-MM. This ELISA kit for research use only, not for therapeutic or diagnostic applications!
Product Categories/Family for CK MM elisa kit	Metabolism; Signal Transduction; Cardiovascular

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/creatine-kinase-mm-isoenzyme-ck-mm/27439>