

General HVD3 (25-Hydroxyvitamin D3) ELISA Kit | HVD3 elisa kit

General HVD3 (25-Hydroxyvitamin D3) ELISA Kit

CATALOG #
AAA275676

REACTIVITY
General


PRODUCT SYNONYMS

HVD3 (25-Hydroxyvitamin D3); General HVD3 (25-Hydroxyvitamin D3) ELISA Kit; 25-HVD3; 25HVD3; Calcifediol; Calcidiol; H-VD3; 25-Hydroxycholecalciferol; 25-Hydroxy Vitamin D3; 25(OH)D; HVD3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	General
Specificity	This assay has high sensitivity and excellent specificity for detection of General HVD3. No significant cross-reactivity or interference between General HVD3 and analogues was observed.
Assay Type	Competitive
Detection Range	3.13-200 ng/mL
Sensitivity	0.91 ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay)? CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.
Inter-assay Precision	Inter-assay Precision (Precision between assays)? CV%<10%. Three samples of known concentration were tested in forty separate assays to assess inter-assay precision.
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HVD3 elisa kit

Principle of the Assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with 25-Hydroxyvitamin D3(HVD3) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to 25-Hydroxyvitamin D3(HVD3). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of 25-Hydroxyvitamin D3(HVD3) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for HVD3 elisa kit

[Reproductive science](#); [Genetic science](#); [Nutrition metabolism](#); [Bone metabolism](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/general-elisa-kits/hvd3-25-hydroxyvitamin-d3/275676>