

Human SLC34A2 (Sodium-dependent phosphate transport protein 2B) ELISA Kit | SLC34A2 elisa kit

Human SLC34A2 (Sodium-dependent phosphate transport protein 2B) ELISA Kit

CATALOG #

AAA27629

REACTIVITY

Human


PRODUCT SYNONYMS

SLC34A2 (Sodium-dependent phosphate transport protein 2B); Human SLC34A2 (Sodium-dependent phosphate transport protein 2B) ELISA Kit; Sodium-dependent phosphate transport protein 2B; Sodium-phosphate transport protein 2B; Na(+)-dependent phosphate cotransporter 2B; NaPi3b; Sodium/phosphate cotransporter 2B; Na(+)/Pi cotransporter 2B; NaPi-2b; Solute carrier family 34 member 2; SLC34A2; SLC34A2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of SLC34A2. No significant cross-reactivity or interference between SLC34A2 and analogues was observed.
Assay Type	Sandwich
Samples	Serum, plasma, tissue homogenates and other biological fluids
Detection Range	31.25-2000pg/ml
Sensitivity	18.75pg/ml
Intra-assay Precision	CV<8%
Inter-assay Precision	CV<10%
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SLC34A2 elisa kit	Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Capture antibody was pre-coated onto 96-well plates. And the biotin conjugated antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and washed with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the target amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of target can be calculated.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	295789158
NCBI GeneID	10568
NCBI Accession #	NP_001171469.1
NCBI GenBank Nucleotide #	NM_001177998.1
UniProt Accession #	O95436; Q8N2K2; Q8WYA9; Q9P0V7; A5PL17
Molecular Weight	75,759 Da
NCBI Official Full Name	sodium-dependent phosphate transport protein 2B isoform b
NCBI Official Synonym Full Names	solute carrier family 34 (type II sodium/phosphate cotransporter), member 2
NCBI Official Symbol	SLC34A2
NCBI Official Synonym Symbols	NPTIIB; NAPI-3B; NAPI-IIB
NCBI Protein Information	sodium-dependent phosphate transport protein 2B; sodium/phosphate cotransporter 2B; type II sodium-dependent phosphate transporter 3b; solute carrier family 34 (sodium phosphate), member 2
UniProt Protein Name	Sodium-dependent phosphate transport protein 2B
UniProt Gene Name	SLC34A2
UniProt Synonym Gene Names	Sodium-phosphate transport protein 2B; Na(+)/Pi cotransporter 2B; NaPi-2b
UniProt Entry Name	NPT2B_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/slc34a2-sodium-dependentphosphatetransportprotein2b/27629>