

General DYKDDDDK-Tag (DYKDDDDK-Tag Protein) ELISA Kit | DYKDDDDK-Tag elisa kit

DYKDDDDK-Tag (DYKDDDDK-Tag Protein) ELISA Kit

CATALOG #
AAA27647

REACTIVITY
General


PRODUCT SYNONYMS

DYKDDDDK-Tag (DYKDDDDK-Tag Protein); DYKDDDDK-Tag (DYKDDDDK-Tag Protein) ELISA Kit; DYKDDDDK-Tag; DYKDDDDK-Tag Protein; DYKDDDDK-Tag elisa kit

PRODUCT SPECIFICATIONS

Reactivity	General
Specificity	This assay has high sensitivity and excellent specificity for detection of FLAG-Tag. No significant cross-reactivity or interference between FLAG-Tag and analogues was observed.
Assay Type	Competitive
Samples	Transfected cell supernatant and lysate, cell lysate, or recombinant protein containing the DYKDDDDK-tag
Detection Range	15.625-1000 ng/ml
Sensitivity	9.375 ng/ml
Preparation and Storage	Store entire kit at 2-8°C for short-term. For longer-term, please store the microplate & standard at -20°C, while the remaining reagents can be stored at 2-8°C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DYKDDDDK-Tag elisa kit

Principle of the Assay: This kit was based on Competitive-ELISA detection method. The microtiter plate provided in this kit has been pre-coated with target. During the reaction, target in the sample or standard competes with a fixed amount of target on the solid phase supporter for sites on the Biotinylated Detection Antibody specific to target. Excess conjugate and unbound sample or standard are washed from the plate, and HRP-Streptavidin (SABC) is added to each microplate well and incubated. Then TMB substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm. The concentration of target in the samples is then determined by comparing the OD of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

