

Cat SDMA (Symmetric dimethylarginine) ELISA Kit | SDMA elisa kit

Cat SDMA (Symmetric dimethylarginine) ELISA Kit

CATALOG #
AAA277019

REACTIVITY
Cat


PRODUCT SYNONYMS

SDMA (Symmetric dimethylarginine); Cat SDMA (Symmetric dimethylarginine) ELISA Kit; Symmetric dimethylarginine; SDMA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Cat
Specificity	This assay has high sensitivity and excellent specificity for detection of Cat SDMA. No significant cross-reactivity or interference between Cat SDMA and analogues was observed.
Assay Type	Competitive
Detection Range	0.16-10 ng/mL
Sensitivity	0.051 ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay)? CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.
Inter-assay Precision	Inter-assay Precision (Precision between assays)? CV%<10%. Three samples of known concentration were tested in forty separate assays to assess inter-assay precision.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SDMA elisa kit

Principle of the Assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with SDMA(Symmetric dimethylarginine) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to SDMA(Symmetric dimethylarginine). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of SDMA(Symmetric dimethylarginine) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for SDMA elisa kit

Signal transduction; Developmental science; Bone metabolism

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.