

## Human CH50 (50% complement hemolysis)ELISA Kit | CH50 elisa kit

Human CH50 (50% complement hemolysis)ELISA Kit

CATALOG #  
**AAA277452**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

CH50 (50% complement hemolysis)ELISA Kit; Human CH50 (50% complement hemolysis)ELISA Kit; CH50; CH50 elisa kit

### PRODUCT SPECIFICATIONS

|                       |                                                                                                                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity            | Human                                                                                                                                                                                 |
| Specificity           | This assay has high sensitivity and excellent specificity for detection of Human CH50. No significant cross-reactivity or interference between Human CH50 and analogues was observed. |
| Assay Type            | Sandwich                                                                                                                                                                              |
| Detection Range       | 2.35-150 U/mL                                                                                                                                                                         |
| Sensitivity           | 1.19 U/mL                                                                                                                                                                             |
| Intra-assay Precision | Intra-assay Precision (Precision within an assay)?CV% < 8%. Three samples of known concentration were tested twenty times on 1 plate to assess intra-assay precision.                 |
| Inter-assay Precision | Inter-assay Precision (precision between assays)?CV% < 10%. Three samples of known concentration were tested in forty separate assays to assess inter-assay precision.                |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for CH50 elisa kit

Principle of the Assay: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human CH50. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Human CH50. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Human CH50, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Human CH50 in the samples is then determined by comparing the OD of the samples to the standard curve.

#### Product Categories/Family for CH50 elisa kit

Immunology; Metabolism; Signal Transduction

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/ch50-50-complement-hemolysis-elisa-kit/277452>