

Rabbit DDT Polyclonal Antibody | anti-DDT antibody

DDT Polyclonal Antibody

CATALOG # AAA296546	HOST Rabbit	REACTIVITY Human, Mouse, Rat	APPLICATIONS Immunohistochemistry
-------------------------------	-----------------------	--	---


PRODUCT SYNONYMS

DDT; Polyclonal Antibody; DDT Polyclonal Antibody; D-dopachrome decarboxylase (EC 4.1.1.84)(D-dopachrome tautomerase)(Phenylpyruvate tautomerase II); anti-DDT antibody


PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Applicable Applications for anti-DDT antibody	IHC (Immunohistochemistry)
Preparation and Storage	Maintain refrigerated at 2-8 degree C for up to 2 weeks. For long term storage store at -20 degree C in small aliquots to prevent freeze-thaw cycles.


IMAGES




IHC (Immunohistochemistry)
Dilution: IHC-p 1:50-200, ELISA 1:10000-20000



IHC (Immunohistochemistry)
Dilution: IHC-p 1:50-200, ELISA 1:10000-20000



IHC (Immunohistochemistry)
Dilution: IHC-p 1:50-200, ELISA 1:10000-20000



IHC (Immunohistochemistry)
Dilution: IHC-p 1:50-200, ELISA 1:10000-20000

ADDITIONAL INFORMATION

Related Product Information for anti-DDT antibody	Tautomerization of D-dopachrome with decarboxylation to give 5, 6-dihydroxyindole (DHI).
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	145386531
NCBI GeneID	1652
NCBI Accession #	NP_001077861.1
NCBI GenBank Nucleotide #	NM_001084392.1
UniProt Accession #	P30046; O00774; O60787; Q13534; B7Z522
Molecular Weight	11,323 Da
NCBI Official Full Name	D-dopachrome decarboxylase
NCBI Official Synonym Full Names	D-dopachrome tautomerase
NCBI Official Symbol	DDT
NCBI Official Synonym Symbols	DDCT
NCBI Protein Information	D-dopachrome decarboxylase; phenylpyruvate tautomerase II
UniProt Protein Name	D-dopachrome decarboxylase
UniProt Gene Name	DDT
UniProt Entry Name	DOPD_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.