

Human Baculoviral IAP repeat-containing protein 2 ELISA Kit | BIRC2 elisa kit

Human Baculoviral IAP repeat-containing protein 2 ELISA Kit

CATALOG #
AAA29751

REACTIVITY
Human


PRODUCT SYNONYMS

Baculoviral IAP repeat-containing protein 2; Human Baculoviral IAP repeat-containing protein 2 ELISA Kit; C-IAP1; IAP homolog B; Inhibitor of apoptosis protein 2; IAP-2; hIAP-2; hIAP2; RING finger protein 48; TNFR2-TRAF-signaling complex protein 2; BIRC2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay recognizes recombinant and natural human BIRC2. No significant cross-reactivity or interference was observed.
Samples	serum, Plasma, tissue homogenates and Cell culture supernates and Other biological fluids
Detection Range	78.0 pg/mL -5000 pg/mL
Preparation and Storage	Store at 4 degree C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BIRC2 elisa kit

Intended Uses: This immunoassay kit allows for the in vitro quantitative determination of human baculoviral iap repeat-containing protein 2,BIRC2 concentrations in serum, Plasma, tissue homogenates and Cell culture supernates and Other biological fluids.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4502141
NCBI GeneID	329
NCBI Accession #	NP_001157.1
NCBI GenBank Nucleotide #	NM_001166.4
UniProt Accession #	Q13490; Q16516; Q4TTG0; B4E026

Molecular Weight	64,098 Da
NCBI Official Full Name	baculoviral IAP repeat-containing protein 2 isoform 1
NCBI Official Synonym Full Names	baculoviral IAP repeat containing 2
NCBI Official Symbol	BIRC2
NCBI Official Synonym Symbols	API1; MIHB; HIAP2; RNF48; cIAP1; Hiap-2; c-IAP1
NCBI Protein Information	baculoviral IAP repeat-containing protein 2
UniProt Protein Name	Baculoviral IAP repeat-containing protein 2
UniProt Gene Name	BIRC2
UniProt Synonym Gene Names	API1; MIHB; RNF48; hIAP-2; hIAP2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/birc2-human-elisa-kits/baculoviral-iap-repeat-containing-protein-2/29751>