

Rabbit anti-Human HICE1 Polyclonal Antibody | anti-HICE1 antibody

HICE1 antibody

CATALOG #	HOST	REACTIVITY	APPLICATIONS
AAA308160	Rabbit	Human	Immunofluorescence, Immunohistochemistry, Western Blot


PRODUCT SYNONYMS

HICE1; Polyclonal Antibody; HICE1 antibody; anti-HICE1 antibody

PRODUCT SPECIFICATIONS


Host	Rabbit
Reactivity	Human
Clonality	Polyclonal
Purity/Purification	Purified by antigen-affinity chromatography.
Form/Format	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Applicable Applications for anti-HICE1 antibody	IF (Immunofluorescence), IHC (Immunohistochemistry), WB (Western Blot)
Target Name	HICE1
Preparation and Storage	Store at -20 degree C for long term preservation (recommended). Store at 4 degree C for short term use.

IMAGES




IF (Immunofluorescence)

Immunofluorescence analysis of human osteosarcoma cell line U2OS, using HICE1 antibody at 1: 50-1: 200 dilution.



IHC (Immunohistochemistry)

A: Immunohistochemical analysis of paraffin-embedded Human lymph tissue, using HICE1 antibody at 1: 50 dilution.(20X) B: Immunohistochemical analysis of paraffin-embedded Human lymph tissue, using HICE1 antibody at 1: 50 dilution.(40X)



Application Data

SampleA: His-Hice1 (2x)B: His-Hice1 (1x)7.5% SDS PAGEPrimary antibody diluted at 1: 500

ADDITIONAL INFORMATION

Related Product Information for anti-HICE1 antibody	HAUS8 is 1 of 8 subunits of the 390-kD human augmin complex, or HAUS complex. The augmin complex was first identified in Drosophila, and its name comes from the Latin verb 'augmentare,' meaning 'to increase.' The augmin complex is a microtubule-binding complex involved in microtubule generation within the mitotic spindle and is vital to mitotic spindle assembly
Product Categories/Family for anti-HICE1 antibody	Total protein Ab

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/polyclonal/hice1-antibody-immunofluorescence-immunohistochemistry-western-blot-elisa/308160>