

Rabbit DNER Polyclonal Antibody | anti-DNER antibody

DNER Rabbit Polyclonal Antibody

CATALOG # AAA313596	HOST Rabbit	REACTIVITY Human, Mouse, Rat	APPLICATIONS Immunofluorescence, Western Blot
-------------------------------	-----------------------	--	---


PRODUCT SYNONYMS

DNER; Polyclonal Antibody; DNER Rabbit Polyclonal Antibody; UNQ26; bet; anti-DNER antibody


PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Isotype	IgG
Purity/Purification	Affinity purification
Form/Format	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Applicable Applications for anti-DNER antibody	IF (Immunofluorescence), WB (Western Blot)
Immunogen	Recombinant fusion protein of human DNER (NP_620711.3).
Preparation and Storage	Store at -20 degree C. Avoid freeze/thaw cycles.

IMAGES


IF (Immunofluorescence)

Immunofluorescence analysis of L929 cells using DNER at dilution of 1:100. Blue: DAPI for nuclear staining.


WB (Western Blot)

Western blot analysis of extracts of various cell lines, using DNER at 1:1000 dilution.

ADDITIONAL INFORMATION

Product Categories/Family for anti-DNER antibody	Primary Antibody
--	------------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	116235485
NCBI GeneID	92737
NCBI Accession #	NP_620711.3
NCBI GenBank Nucleotide #	NM_139072.3
UniProt Accession #	Q8NFT8; Q53R88; Q53TP7; Q53TQ5; Q8IYT0; Q8TB42; Q9NTF1; Q9UDM2; A6NP39
Molecular Weight	~120kD (SDS-PAGE)
NCBI Official Full Name	delta and Notch-like epidermal growth factor-related receptor
NCBI Official Synonym Full Names	delta/notch-like EGF repeat containing
NCBI Official Symbol	DNER
NCBI Official Synonym Symbols	bet; UNQ26
NCBI Protein Information	delta and Notch-like epidermal growth factor-related receptor; H_NH0150O02.1; WUGSC:H_NH0150O02.1; delta-notch-like EGF repeat-containing transmembrane
UniProt Protein Name	Delta and Notch-like epidermal growth factor-related receptor
UniProt Gene Name	DNER
UniProt Synonym Gene Names	BET
UniProt Entry Name	DNER_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/polyclonal/dner-antibody-immunofluorescence-western-blot-elisa/313596>