

Human, Mouse TIP60 ELISA Kit | TIP60 elisa kit

Phospho-TIP60 (Ser86) Colorimetric Cell-Based ELISA Kit

CATALOG # AAA315577	REACTIVITY Human, Mouse
-------------------------------	-----------------------------------

PRODUCT SYNONYMS

TIP60; Phospho-TIP60 (Ser86) Colorimetric Cell-Based ELISA Kit; 60 kDa Tat interactive protein; cPLA(2)-interacting protein; histone acetyltransferase HTATIP; HIV-1 Tat interactive protein; HTATIP; Tat-interactive 60 kDa protein alpha isoform; Tat-interactive protein-60; TI60; TIP60 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human, Mouse
Specificity	The Anti-TIP60 (Phospho-Ser86) Antibody is a rabbit polyclonal antibody. It was tested on Western Blots for specificity. The data in Figure 2 shows that a single protein band was detected. This protein band can be blocked by the synthesized immunogen peptide.
Form/Format	Two 96-Well Plates
Sequence Length	494
Samples	Researchers plate their cell line of choice
Assay Type	Phospho Cell-Based ELISA
Sensitivity	>=5000 cells/well
Reactivity Note	Human: S86, Mouse: S86
Detection Method	Colorimetric 450 nm
Preparation and Storage	Store at 4 degree C.

ADDITIONAL INFORMATION

Related Product Information for TIP60 elisa kit	Background: The Colorimetric Cell-Based ELISA Kit allows for the detection of various target proteins and the effects that certain stimulation conditions have on target protein expression in different cell lines. Qualitative determination of target protein concentration is achieved by an indirect ELISA format. In essence, the target protein is captured by target-specific primary (1 degree) antibodies while the HRP-conjugated secondary (2 degree) antibodies bind the Fc region of the 1 degree antibody. Through this binding, the HRP enzyme conjugated to the 2 degree antibody can catalyze a colorimetric reaction upon substrate addition. Due to the qualitative nature of the CellBased ELISA, multiple normalization methods are described: 1) a monoclonal antibody specific for human GAPDH is included to serve as an internal positive control in normalizing the target absorbance values. 2) Following the colorimetric measurement of HRP activity via substrate addition, the Crystal Violet whole-cell staining method is used to determine cell density. After staining, the results can be analyzed by normalizing the absorbance values to cell amounts, by which the plating difference can be adjusted. 3) If a phosphorylated target is being detected, an antibody against the non- phosphorylated counterpart will be provided for normalization purposes. The absorbance values obtained for the nonphosphorylated target can be used to normalize the absorbance values for the phosphorylated target.
Product Categories/Family for TIP60 elisa kit	Cell Based ELISA Kit

NCBI AND UNIPROT INFORMATION

NCBI GI #	332205931
-----------	-----------

NCBI GeneID	10524
NCBI Accession #	NP_001193762.1
NCBI GenBank Nucleotide #	NM_001206833.1
UniProt Accession #	Q92993; O95624; Q13430; Q17RW5; Q561W3; Q6GSE8; Q9BWK7; B4E3C7; C9JL99
Molecular Weight	58582
NCBI Official Full Name	histone acetyltransferase KAT5 isoform 4
NCBI Official Synonym Full Names	lysine acetyltransferase 5
NCBI Official Symbol	KAT5
NCBI Official Synonym Symbols	TIP; ESA1; PLIP; TIP60; cPLA2; HTATIP; ZC2HC5; HTATIP1
NCBI Protein Information	histone acetyltransferase KAT5
UniProt Protein Name	Histone acetyltransferase KAT5
UniProt Gene Name	KAT5
UniProt Synonym Gene Names	HTATIP; TIP60; Tip60; HIV-1 Tat interactive protein

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-mouse-elisa-kits/tip60/315577>