

Human RANTES ELISA Kit | RANTES elisa kit

RANTES (Human) OmniKine ELISA Kit

CATALOG #
AAA315649

REACTIVITY
Human


PRODUCT SYNONYMS

RANTES; RANTES (Human) OmniKine ELISA Kit; C-C motif chemokine 5; EoCP; Eosinophil chemotactic cytokine; SIS-delta; Small-inducible cytokine A5; T cell-specific protein P228; T-cell-specific protein RANTES; D17S136E; SCYA5; RANTES elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The OmniKine Human RANTES ELISA is capable of recognizing both recombinant and naturally produced Human RANTES proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference. Human: Eotaxin, LEC, MCP-2, MCP-3, MCP-4, MIP-1alpha, MIP-1beta, MIP-5 Murine: MCP-2, MCP-3, MIP-1alpha, RANTES Rat: RANTES
Form/Format	12 x 8-Well Microstrips
Sequence Length	91
Samples	Plasma, serum, lysate
Assay Type	Sandwich ELISA
Detection Range	16-2000 pg/ml
Sensitivity	16-2000 pg/ml
Detection Method	Colorimetric 450 nm
Subtype	None
Preparation and Storage	Store entire kit at 4 degree C until use.

IMAGES



Application Data

The OmniKine Human RANTES ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human RANTES proteins within the range of 16-2000 pg/ml.

NCBI AND UNIPROT INFORMATION

NCBI GI #	22538814
NCBI GeneID	6352
NCBI Accession #	NP_002976.2
NCBI GenBank Nucleotide #	NM_002985.2
UniProt Accession #	P13501; O43646; Q0QVW8; Q4ZGJ1; Q9NYA2; Q9UBG2; Q9UC99
Molecular Weight	9,990 Da
NCBI Official Full Name	C-C motif chemokine 5 isoform 1
NCBI Official Synonym Full Names	C-C motif chemokine ligand 5
NCBI Official Symbol	CCL5
NCBI Official Synonym Symbols	SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E; SIS-delta
NCBI Protein Information	C-C motif chemokine 5
UniProt Protein Name	C-C motif chemokine 5
UniProt Gene Name	CCL5
UniProt Synonym Gene Names	D17S136E; SCYA5; TCP228

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/rantes/315649>