

## Human TACI ELISA Kit | TACI elisa kit

TACI (Human) OmniKine ELISA Kit

CATALOG #  
**AAA315654**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

TACI; TACI (Human) OmniKine ELISA Kit; Tumor necrosis factor receptor superfamily member 13B; Transmembrane activator and CAML interactor; CD267; TACI elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Human                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Specificity             | The OmniKine Human TACI ELISA is capable of recognizing both recombinant and naturally produced Human TACI proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference.<br>Human: 4-1BBL, 4-1BBR, AITRL, BAFF, BAFF-R, BCMA, CD-40L, Fas-L, Fas-R, OPG, sRANK-L, sRANK-R, TNF-alpha, TNF-beta, sTNF-R1, sTNF-R2, sTRAIL, sTRAIL-R1, sTRAIL-R2, TWEAK, TWEAK-R |
| Form/Format             | 12 x 8-Well Microstrips                                                                                                                                                                                                                                                                                                                                                                                                        |
| Sequence Length         | 293                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Samples                 | Plasma, serum, lysate                                                                                                                                                                                                                                                                                                                                                                                                          |
| Assay Type              | Sandwich ELISA                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Detection Range         | 16-2000 pg/ml                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Sensitivity             | 16-2000 pg/ml                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Detection Method        | Colorimetric 450 nm                                                                                                                                                                                                                                                                                                                                                                                                            |
| Subtype                 | None                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Preparation and Storage | Store entire kit at 4 degree C until use.                                                                                                                                                                                                                                                                                                                                                                                      |

### IMAGES



#### Application Data

The OmniKine Human TACI ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human TACI proteins within the range of 16-2000 pg/ml.

NCBI AND UNIPROT INFORMATION

|                                  |                                                       |
|----------------------------------|-------------------------------------------------------|
| NCBI GI #                        | 6912694                                               |
| NCBI GeneID                      | 23495                                                 |
| NCBI Accession #                 | NP_036584.1                                           |
| NCBI GenBank Nucleotide #        | NM_012452.2                                           |
| UniProt Accession #              | O14836; Q32LX4; Q7Z6F5; B2R8B0; B7Z6V8                |
| Molecular Weight                 | 19,865 Da                                             |
| NCBI Official Full Name          | tumor necrosis factor receptor superfamily member 13B |
| NCBI Official Synonym Full Names | TNF receptor superfamily member 13B                   |
| NCBI Official Symbol             | TNFRSF13B                                             |
| NCBI Official Synonym Symbols    | CVID; RYZN; TACI; CD267; CVID2; IGAD2; TNFRSF14B      |
| NCBI Protein Information         | tumor necrosis factor receptor superfamily member 13B |
| UniProt Protein Name             | Tumor necrosis factor receptor superfamily member 13B |
| UniProt Gene Name                | <b>TNFRSF13B</b>                                      |
| UniProt Synonym Gene Names       | TACI                                                  |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/taci/315654>