

Human TWEAK ELISA Kit | TWEAK elisa kit

TWEAK (Human) OmniKine ELISA Kit

CATALOG #
AAA315661

REACTIVITY
Human


PRODUCT SYNONYMS

TWEAK; TWEAK (Human) OmniKine ELISA Kit; Tumor necrosis factor ligand superfamily member 12; APO3 ligand; TNF-related weak inducer of apoptosis; APO3L; DR3LG; TWEAK elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The OmniKine Human TWEAK ELISA is capable of recognizing both recombinant and naturally produced Human TWEAK proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference. Human: 4-1BBL, 4-1BBR, BAFF, CD40L, sFasL, OPG, sRANKL, sRANKR, TNFalpha, TNFbeta, sTNF-RI, sTNF-RII, sTRAIL, sTRAIL-RI, sTRAIL-RII, TWEAKR Murine: sCD40L, sRANKL, TNFalpha Rat: TNFalpha
Form/Format	12 x 8-Well Microstrips
Sequence Length	249
Samples	Plasma, serum, lysate
Assay Type	Sandwich ELISA
Detection Range	16-2000 pg/ml
Sensitivity	16-2000 pg/ml
Detection Method	Colorimetric 450 nm
Subtype	None
Preparation and Storage	Store entire kit at 4 degree C until use.

IMAGES



Application Data

The OmniKine Human TWEAK ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human TWEAK proteins within the range of 16-2000 pg/ml.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507597
NCBI GeneID	8742
NCBI Accession #	NP_003800.1
NCBI GenBank Nucleotide #	NM_003809.2
UniProt Accession #	O43508; Q8IZK7; Q8WUZ7
Molecular Weight	36,589 Da
NCBI Official Full Name	tumor necrosis factor ligand superfamily member 12 proprotein
NCBI Official Synonym Full Names	TNF superfamily member 12
NCBI Official Symbol	TNFSF12
NCBI Official Synonym Symbols	APO3L; DR3LG; TWEAK; TNLG4A
NCBI Protein Information	tumor necrosis factor ligand superfamily member 12
UniProt Protein Name	Tumor necrosis factor ligand superfamily member 12
UniProt Gene Name	TNFSF12
UniProt Synonym Gene Names	APO3L; DR3LG; TWEAK

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tweak/315661>