

Human FasR ELISA Kit | FasR elisa kit

FasR (Human) OmniKine ELISA Kit

CATALOG #
AAA315709

REACTIVITY
Human


PRODUCT SYNONYMS

FasR; FasR (Human) OmniKine ELISA Kit; Tumor necrosis factor receptor superfamily member 6; Apo-1 antigen; Apoptosis-mediating surface antigen FAS; FASLG receptor; CD95; APT1; FAS1; TNFRSF6; FasR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The OmniKine Human Fas Receptor ELISA is capable of recognizing both recombinant and naturally produced Human Fas Receptor proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference. Human: Fas-Ligand, RANTES, TGF-alpha, TNF-alpha Murine: FasR, Fas-Ligand, TNF-R1, TNF-R2
Form/Format	12 x 8-Well Microstrips
Sequence Length	335
Samples	Plasma, serum, lysate
Assay Type	Sandwich ELISA
Detection Range	63-4000 pg/ml
Sensitivity	63-4000 pg/ml
Detection Method	Colorimetric 450 nm
Subtype	None
Preparation and Storage	Store entire kit at 4 degree C until use.

IMAGES



Application Data

The OmniKine Human FasR ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human FasR proteins within the range of 63-4000 pg/ml.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507583
NCBI GeneID	355
NCBI Accession #	NP_000034.1
NCBI GenBank Nucleotide #	NM_000043.5
UniProt Accession #	P25445; Q14292; Q14293; Q14294; Q14295; Q16652; Q5T9P1; Q5T9P2; Q5T9P3; Q6SSE9; A9UJX4; B6VNV4
Molecular Weight	24,781 Da
NCBI Official Full Name	tumor necrosis factor receptor superfamily member 6 isoform 1
NCBI Official Synonym Full Names	Fas cell surface death receptor
NCBI Official Symbol	FAS
NCBI Official Synonym Symbols	APT1; CD95; FAS1; APO-1; FASTM; ALPS1A; TNFRSF6
NCBI Protein Information	tumor necrosis factor receptor superfamily member 6
UniProt Protein Name	Tumor necrosis factor receptor superfamily member 6
UniProt Gene Name	FAS
UniProt Synonym Gene Names	APT1; FAS1; TNFRSF6

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/fasr/315709>