

Human VEGF-R3 ELISA Kit | VEGF-R3 elisa kit

VEGF-R3 (Human) OmniKine ELISA Kit

CATALOG #
AAA315749

REACTIVITY
Human


PRODUCT SYNONYMS

VEGF-R3; VEGF-R3 (Human) OmniKine ELISA Kit; Vascular endothelial growth factor receptor 3; Fms-like tyrosine kinase 4; Tyrosine-protein kinase receptor FLT4; VEGF-R3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The OmniKine Human VEGF-R3 ELISA is capable of recognizing both recombinant and naturally produced Human VEGF-R3 proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference. Human: VEGF121, VEGF165, MCP-2, VEGF-R1/Fc Chimera, VEGF-R2/Fc Chimera
Form/Format	12 x 8-Well Microstrips
Sequence Length	1298
Samples	Plasma, serum, lysate
Assay Type	Sandwich ELISA
Detection Range	157-10000 pg/ml
Sensitivity	157-10000 pg/ml
Detection Method	Colorimetric 450 nm
Subtype	None
Preparation and Storage	Store entire kit at 4 degree C until use.

IMAGES



Application Data

The OmniKine Human VEGF-R3 ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human VEGF-R3 proteins within the range of 157-10000 pg/ml.

NCBI AND UNIPROT INFORMATION

NCBI GI #	103472027
NCBI GeneID	2324
NCBI Accession #	NP_002011.2
NCBI GenBank Nucleotide #	NM_002020.4
UniProt Accession #	P35916; Q16067; Q86W07; Q86W08; A8K6L4; B5A926
Molecular Weight	93,174 Da
NCBI Official Full Name	vascular endothelial growth factor receptor 3 isoform 2
NCBI Official Synonym Full Names	fms related tyrosine kinase 4
NCBI Official Symbol	FLT4
NCBI Official Synonym Symbols	PCL; FLT-4; FLT41; LMPH1A; VEGFR3; VEGFR-3
NCBI Protein Information	vascular endothelial growth factor receptor 3
UniProt Protein Name	Vascular endothelial growth factor receptor 3
UniProt Gene Name	FLT4
UniProt Synonym Gene Names	VEGFR3; VEGFR-3; FLT-4

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/vegfr3/315749>