

Human uPAR ELISA Kit | uPAR elisa kit

uPAR (Human) OmniKine ELISA Kit

CATALOG #
AAA315842

REACTIVITY
Human


PRODUCT SYNONYMS

uPAR; uPAR (Human) OmniKine ELISA Kit; Urokinase plasminogen activator surface receptor; U-PAR; Monocyte activation antigen Mo3; CD87; MO3; UPAR; uPAR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The OmniKine Human uPAR ELISA is capable of recognizing both recombinant and naturally produced Human uPAR proteins. The antigens listed below were tested and interfere at a concentration greater than 50 ng/ml but did not exhibit significant cross reactivity. Human: uPA The antigens listed below were tested and interfere at a concentration greater than 1.56 ng/ml but did not exhibit significant cross reactivity. Murine: uPAR
Form/Format	12 x 8-Well Microstrips
Sequence Length	281
Samples	Plasma, serum, lysate
Assay Type	Sandwich ELISA
Detection Range	32-2000 pg/ml
Sensitivity	32-2000 pg/ml
Detection Method	Colorimetric 450 nm
Subtype	None
Preparation and Storage	Store entire kit at 4 degree C until use.

IMAGES



Application Data

The OmniKine Human uPAR ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human uPAR proteins within the range of 32-2000 pg/ml.

NCBI AND UNIPROT INFORMATION

NCBI GI #	53829379
NCBI GeneID	5329
NCBI Accession #	NP_001005376.1
NCBI GenBank Nucleotide #	NM_001005376.2
UniProt Accession #	Q03405; Q12876; Q15845; Q16887; Q6IB52; Q9BWT0; Q9NYC8; Q9UD69; Q9UEA6; Q9UM92; Q9UMV0; A8K409
Molecular Weight	32,016 Da
NCBI Official Full Name	urokinase plasminogen activator surface receptor isoform 2
NCBI Official Synonym Full Names	plasminogen activator, urokinase receptor
NCBI Official Symbol	PLAUR
NCBI Official Synonym Symbols	CD87; UPAR; URKR; U-PAR
NCBI Protein Information	urokinase plasminogen activator surface receptor
UniProt Protein Name	Urokinase plasminogen activator surface receptor
UniProt Gene Name	PLAUR
UniProt Synonym Gene Names	MO3; UPAR; U-PAR; uPAR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/upar/315842>