

Human BD-4 CLIA Kit | BD-4 clia kit

BD-4 (Human) LumiAb ELISA Kit

CATALOG #
AAA318520

REACTIVITY
Human


PRODUCT SYNONYMS

BD-4; BD-4 (Human) LumiAb ELISA Kit; Human BD-4; HBD-4; DEFB-4; HBD4; DEFB104B; Beta-defensin 4; BD-4 clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The LumiAb Human BD-3 ELISA is capable of recognizing both recombinant and naturally produced Human BD-3 proteins. The antigens listed below were tested at 50 ng/ml and did not exhibit significant cross reactivity or interference.
Sequence Length	72
Assay Type	Quantitative Sandwich
Sensitivity	15 pg/ml.
Preparation and Storage	Store at 4 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BD-4 clia kit	Principle of the Assay: The LumiAb Human BD-3 Chemiluminescent ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Human BD-3 concentrations within any experimental sample including cell lysates, serum and plasma. This particular immunoassay utilizes the quantitative technique of a "Sandwich" Enzyme-Linked Immunosorbent Assay (ELISA) where the target protein (antigen) is bound in a "sandwich" format by the primary capture antibodies coated to each well-bottom and the secondary detection antibodies added subsequently by the investigator. The capture antibodies coated to the bottom of each well are specific for a particular epitope on Human BD-3 while the user-added detection antibodies bind to epitopes on the captured target protein. Amid each step of the procedure, a series of wash steps must be performed to ensure the elimination of non-specific binding between proteins to other proteins or to the solid phase. After incubation and "sandwiching" of the target antigen, a peroxidase enzyme is conjugated to the constant heavy chain of the secondary antibody (either covalently or via Avidin/Streptavidin-Biotin interactions), allowing for a sensitive luminescent reaction to ensue upon substrate addition. When the Peroxide Enhancer solution is added, the reaction catalyzed by peroxidase yields light that is representative of the antigen concentration. After a brief incubation, the microplate can be read with a luminometer, allowing for generation of a standard curve and subsequent determination of protein concentration!!Background/Introduction: Human BD-3, also known as Beta-Defensin 3, is a 67 amino acid cytokine protein encoded by the DEFB103A gene located at locus 8p23.1 on chromosome 8. After signal sequence processing and maturation, the cytokine becomes 45 amino acid residues long (from 23 to 67). Human BD-3 is known to exhibit antimicrobial activity against Gram-positive bacteria, S.aureus and S.pyogenes, Gram-negative bacteria P.aeruginosa and E.coli and the yeast C.albicans. Moreover, the cytokine kills multiresistant S.aureus and vancomycin-resistant E.faecium. This particular BD protein is expressed highly in skin and tonsils, and to a lesser extent in trachea, uterus, kidney, thymus, adenoid, pharynx and tongue. It has much lower expression in salivary glands, bone marrow, colon, stomach, polyps and larynx along with no expression in the small intestine. Furthermore, induction of BD-3 is achieved through bacterial infection or detection of lipopolysaccharides (LPS) and by IFN-gamma.
Product Categories/Family for BD-4 clia kit	Luminometer Sandwich ELISA

NCBI AND UNIPROT INFORMATION

NCBI GI #	103485137
NCBI GeneID	503618
NCBI Accession #	NP_001035792.1
NCBI GenBank Nucleotide #	NM_001040702.1
UniProt Accession #	Q8WTQ1
NCBI Official Full Name	beta-defensin 104
NCBI Official Synonym Full Names	defensin beta 104B
NCBI Official Symbol	DEFB104B
NCBI Official Synonym Symbols	BD-4; hBD-4; DEFB-4
NCBI Protein Information	beta-defensin 104
UniProt Protein Name	Beta-defensin 104
UniProt Gene Name	DEFB104A
UniProt Synonym Gene Names	DEFB104; DEFB4; BD-4; DEFB-4; hBD-4
UniProt Entry Name	D104A_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/bd-4/318520>