

Bovine Bcl-2 ELISA Kit | Bcl-2 elisa kit

Bcl-2 (Phospho-Ser70) Colorimetric Cell-Based ELISA Kit

CATALOG #
AAA319291

REACTIVITY
Bovine


PRODUCT SYNONYMS

Bcl-2; Bcl-2 (Phospho-Ser70) Colorimetric Cell-Based ELISA Kit; Bcl-2 (Phospho-Ser70) Cell-Based ELISA Kit; B-cell differentiation factor; B-cell stimulatory factor 2; BSF2; BSF-2; CDF; CTL differentiation factor ; HSF; hybridoma growth factor; IFNB2; IFN-beta-2; IL6; IL-6; Interferon beta-2; interleukin 6 (interferon; beta 2); interleukin BSF-2; interleukin-6; MGI-2A; Bcl-2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Form/Format	Sandwich ELISA
Samples	Plasma, Serum, Cell Culture Supernatant
Detection Range	16-1,000 pg/mL
Detection Method	Colorimetric 450 nm
Preparation and Storage	4 degree C/3 Months

IMAGES



WB (Western Blot)

Western Blot analysis of JK cells using Phospho-Bcl-2 (S70) Polyclonal Antibody

ADDITIONAL INFORMATION

Related Product Information for Bcl-2 elisa kit

The Bcl-2 (Phospho-Ser70) Colorimetric Cell-Based ELISA Kit is a convenient, lysate-free, high throughput and sensitive assay kit that can detect Bcl-2 (Phospho-Ser70) protein expression profile in cells. The kit can be used for measuring the relative amounts of Bcl-2 (Phospho-Ser70) in cultured cells as well as screening for the effects that various treatments, inhibitors (ie siRNA or chemicals), or activators have on Bcl-2 (Phospho-Ser70). This kit contains Crystal Violet Staining normalization protocol. It provides a means to normalize Absorbance 450 nm of primary antibodies to cell density of the cells in the same well. This provides an accurate account for stimulations and cell densities during the assay. By comparing the proportional values obtained from the experimental wells to those of the control wells (e.g. untreated cells), it is possible to determine the extent to which the stimulant(s) have affected the expression levels.

NCBI AND UNIPROT INFORMATION

NCBI GeneID	280826
-------------	--------

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/bovine-ovine-elisa-kits/bcl-2/319291>