

Rabbit TAK1 Polyclonal Antibody | anti-TAK1 antibody

Phospho-TAK1 (Ser439) Antibody

CATALOG #	HOST	REACTIVITY	APPLICATIONS
AAA321401	Rabbit	Human, Mouse, Rat	Immunofluorescence, Immunohistochemistry, Western Blot


PRODUCT SYNONYMS

TAK1; Polyclonal Antibody; Phospho-TAK1 (Ser439) Antibody; M3K7_HUMAN; MAP3K 7; Map3k7; MEKK7; Mitogen activated protein kinase kinase kinase 7; Mitogen-activated protein kinase kinase kinase 7; TGF beta activated kinase 1; TGF-beta-activated kinase 1; TGF1a; Transforming growth factor beta activated kinase 1; Transforming growth factor-beta-activated kinase 1; anti-TAK1 antibody

PRODUCT SPECIFICATIONS


Host	Rabbit
Reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Isotype	IgG
Specificity	Phospho-TAK1 (Ser439) Antibody detects endogenous levels of TAK1 only when phosphorylated at Ser439.
Purity/Purification	Peptide affinity purification.
Form/Format	Liquid Phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1mg/ml (varies by lot)
Sequence Length	579
Applicable Applications for anti-TAK1 antibody	IF (Immunofluorescence), IHC (Immunohistochemistry), WB (Western Blot)
Immunogen	A synthesized peptide derived from human TAK1 around the phosphorylation site of Ser439.
Preparation and Storage	Store at -20 degree C. Stable for 12 months from date of receipt.

IMAGES



IF (Immunofluorescence)

AAA321401 staining HepG2 cells (30min of 4uM Forskolin treatment) by IF/ICC. The samples were fixed with PFA and permeabilized in 0.1% Triton X-100, then blocked in 10% serum for 45 minutes at 25°C. Samples were then incubated with primary Ab(AAA321401) and mouse anti-beta tubulin Ab for 1 hour at 37°C. An AlexaFluor594 conjugated goat anti-rabbit IgG(H+L) Ab(Red) and an AlexaFluor488 conjugated goat anti-mouse IgG(H+L) Ab(Green) were used as the secondary antibody. The nuclear counter stain is DAPI(blue).



IHC (Immunohistochemistry)

AAA321401 at 1/100 staining Human thyroid cancer tissue by IHC-P. The sample was formaldehyde fixed and a heat mediated antigen retrieval step in citrate buffer was performed. The sample was then blocked and incubated with the antibody for 1.5 hours at 22°C. An HRP conjugated goat anti-rabbit antibody was used as the secondary.

ADDITIONAL INFORMATION

Related Product Information for anti-TAK1 antibody	Serine/threonine kinase which acts as an essential component of the MAP kinase signal transduction pathway. Plays an important role in the cascades of cellular responses evoked by changes in the environment. Mediates signal transduction of TRAF6, various cytokines including interleukin-1 (IL-1), transforming growth factor-beta (TGFB), TGFB-related factors like BMP2 and BMP4, toll-like receptors (TLR), tumor necrosis factor receptor CD40 and B-cell receptor (BCR). Ceramides are also able to activate MAP3K7/TAK1. Once activated, acts as an upstream activator of the MKK/JNK signal transduction cascade and the p38 MAPK signal transduction cascade through the phosphorylation and activation of several MAP kinase kinases like MAP2K1/MEK1, MAP2K3/MKK3, MAP2K6/MKK6 and MAP2K7/MKK7.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507361
NCBI GeneID	6885
NCBI Accession #	NP_003179.1
NCBI GenBank Nucleotide #	NM_003188.3
UniProt Accession #	O43318; O43317; O43319; Q5TDN2; Q5TDN3; Q5TDT7; Q9NTR3; Q9NZ70; B2RE27; E1P523
Molecular Weight	82 kDa
NCBI Official Full Name	mitogen-activated protein kinase kinase kinase 7 isoform A
NCBI Official Synonym Full Names	mitogen-activated protein kinase kinase kinase 7
NCBI Official Symbol	MAP3K7
NCBI Official Synonym Symbols	CSCF; FMD2; TAK1; MEKK7; TGF1a
NCBI Protein Information	mitogen-activated protein kinase kinase kinase 7
UniProt Protein Name	Mitogen-activated protein kinase kinase kinase 7
UniProt Gene Name	MAP3K7
UniProt Synonym Gene Names	TGF-beta-activated kinase 1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.