

Rat Thyroid-Peroxidase (TPO) ELISA Kit | TPO elisa kit

Rat Thyroid-Peroxidase (TPO) ELISA Kit

CATALOG #
AAA330172

REACTIVITY
Rat


PRODUCT SYNONYMS

Thyroid-Peroxidase (TPO); Rat Thyroid-Peroxidase (TPO) ELISA Kit; TPO; MSA; TDH2A; TPX; OTTHUMP00000115536; thyroid microsomal antigen; thyroperoxidase; TPO elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Specificity	Rat Thyroid-Peroxidase (TPO) ELISA Kit has high sensitivity and excellent specificity for detection of Rat TPO. No significant cross-reactivity or interference between Rat TPO and analogues was observed.
Sequence Length	914
Assay Type	Sandwich ELISA (Quantitative)
Samples	Cell Culture Supernatants, Serum, Plasma, Other Biological Fluids
Detection Method	Colorimetric
Usage Notes	* Do not mix components from different kit lots or use reagents beyond the kit expiration date. * Allow all reagents to warm to room temperature for at least 30 minutes before opening. * Pre-rinse the pipet tip with reagent, use fresh pipet tips for each sample, standard and reagent to avoid contamination. * Unused wells must be kept desiccated at 4 degree C in the sealed bag provided. * Mix Thoroughly is very important for the result. It is recommended using low frequency oscillator or slight hand shaking every 10 minutes. * It is recommended that all samples and standards be assayed in duplicate or triplicate.
Preparation and Storage	Store at 2-8 degree C.

IMAGES



Standard Curve (Sample)

Fig.1. Rat Thyroid-Peroxidase (TPO) Standard Curve.

ADDITIONAL INFORMATION

Related Product Information for TPO elisa kit	This Rat Thyroid-Peroxidase (TPO) ELISA Kit employs a two-site sandwich ELISA to quantitate TPO in samples. An antibody specific for TPO has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TPO present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TPO is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TPO bound in the initial step. The color development is stopped and the intensity of the color is measured.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	129832
NCBI GeneID	54314
UniProt Accession #	P14650
NCBI Official Full Name	Thyroid peroxidase
NCBI Official Synonym Full Names	thyroid peroxidase
NCBI Official Symbol	Tpo
NCBI Protein Information	thyroid peroxidase
UniProt Protein Name	Thyroid peroxidase
UniProt Gene Name	Tpo
UniProt Synonym Gene Names	TPO

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/thyroid-peroxidase-tpo/330172>