

Human Tumor Necrosis Factor-Related Apoptosis-Inducing Ligand Receptor 1 (TRAIL-R1) ELISA Kit | TRAIL-R1 elisa kit

Human Tumor Necrosis Factor-Related Apoptosis-Inducing Ligand Receptor 1 (TRAIL-R1) ELISA Kit

CATALOG #
AAA330447

REACTIVITY
Human


PRODUCT SYNONYMS

Tumor Necrosis Factor-Related Apoptosis-Inducing Ligand Receptor 1 (TRAIL-R1); Human Tumor Necrosis Factor-Related Apoptosis-Inducing Ligand Receptor 1 (TRAIL-R1) ELISA Kit; TNFRSF10A; APO2; CD261; DR4; MGC9365; TRAILR-1; TRAILR1; TNF-related apoptosis inducing ligand receptor 1; cytotoxic TRAIL receptor; death receptor 4; TRAIL-R1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	Human Tumor necrosis factor-related apoptosis-inducing ligand receptor 1 (TRAIL-R1) ELISA Kit has high sensitivity and excellent specificity for detection of Human TNFRSF10A. No significant cross-reactivity or interference between Human TNFRSF10A and analogues was observed.
Sequence Length	468
Assay Type	Sandwich ELISA (Quantitative)
Samples	Cell Culture Supernatants, Serum, Plasma, Other Biological Fluids
Detection Method	Colorimetric
Usage Notes	* Do not mix components from different kit lots or use reagents beyond the kit expiration date. * Allow all reagents to warm to room temperature for at least 30 minutes before opening. * Pre-rinse the pipet tip with reagent, use fresh pipet tips for each sample, standard and reagent to avoid contamination. * Unused wells must be kept desiccated at 4 degree C in the sealed bag provided. * Mix Thoroughly is very important for the result. It is recommended using low frequency oscillator or slight hand shaking every 10 minutes. * It is recommended that all samples and standards be assayed in duplicate or triplicate.
Preparation and Storage	Store at 2-8 degree C.

IMAGES



Standard Curve (Sample)

Fig.1. Human Tumor necrosis factor-related apoptosis-inducing ligand receptor 1 (TRAIL-R1) Standard Curve.

ADDITIONAL INFORMATION

Related Product Information for TRAIL-R1 elisa kit	Description: This Human Tumor necrosis factor-related apoptosis-inducing ligand receptor 1 (TRAIL-R1) ELISA Kit employs a two-site sandwich ELISA to quantitate TNFRSF10A in samples. An antibody specific for TNFRSF10A has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNFRSF10A present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNFRSF10A is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNFRSF10A bound in the initial step. The color development is stopped and the intensity of the color is measured. Background: Tumor necrosis factor-related apoptosis-inducing ligand receptor 1 (TRAIL-R1) is a member of the TNF-receptor superfamily. TRAIL-R1 is activated by tumor necrosis factor-related apoptosis inducing ligand (TNFSF10/TRAIL), and thus transduces cell death signal and induces cell apoptosis. Studies with FADD-deficient mice suggested that FADD, a death domain containing adaptor protein, is required for the apoptosis mediated by TRAIL-R1.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	259906438
NCBI GeneID	8797
NCBI Accession #	NP_003835.3
NCBI GenBank Nucleotide #	NP_003835.3
UniProt Accession #	O00220
NCBI Official Full Name	tumor necrosis factor receptor superfamily member 10A
NCBI Official Synonym Full Names	TNF receptor superfamily member 10a
NCBI Official Symbol	TNFRSF10A
NCBI Official Synonym Symbols	DR4; APO2; CD261; TRAILR1; TRAILR-1
NCBI Protein Information	tumor necrosis factor receptor superfamily member 10A
UniProt Protein Name	Tumor necrosis factor receptor superfamily member 10A
UniProt Gene Name	TNFRSF10A
UniProt Synonym Gene Names	APO2; DR4; TRAILR1; TRAIL receptor 1; TRAIL-R1
UniProt Entry Name	TR10A_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.