

Mouse Small Ubiquitin-Related Modifier 1 (SUMO1) ELISA Kit | SUMO1 elisa kit

Mouse Small Ubiquitin-Related Modifier 1 (SUMO1) ELISA Kit

CATALOG #
AAA332292

REACTIVITY
Mouse


PRODUCT SYNONYMS

Small Ubiquitin-Related Modifier 1 (SUMO1); Mouse Small Ubiquitin-Related Modifier 1 (SUMO1) ELISA Kit; SUMO1; DAP-1; GMP1; OFC10; PIC1; SENP2; SMT3; SMT3C; SMT3H3; SUMO-1; UBL1; GAP modifying protein 1; SMT3 suppressor of mif two 3 homolog 1; ubiquitin-like 1 (sentrin); SUMO1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	Mouse Small ubiquitin-related modifier 1 (SUMO1) ELISA Kit has high sensitivity and excellent specificity for detection of Mouse SUMO1. No significant cross-reactivity or interference between Mouse SUMO1 and analogues was observed.
Sequence Length	101
Assay Type	Sandwich ELISA (Quantitative)
Samples	Cell Culture Supernatants, Serum, Plasma, Other Biological Fluids
Detection Method	Colorimetric
Usage Notes	* Do not mix components from different kit lots or use reagents beyond the kit expiration date. * Allow all reagents to warm to room temperature for at least 30 minutes before opening. * Pre-rinse the pipet tip with reagent, use fresh pipet tips for each sample, standard and reagent to avoid contamination. * Unused wells must be kept desiccated at 4 degree C in the sealed bag provided. * Mix Thoroughly is very important for the result. It is recommended using low frequency oscillator or slight hand shaking every 10 minutes. * It is recommended that all samples and standards be assayed in duplicate or triplicate.
Preparation and Storage	Store at 2-8 degree C.

IMAGES



Standard Curve (Sample)

Fig.1. Mouse Small ubiquitin-related modifier 1 (SUMO1) Standard Curve.

ADDITIONAL INFORMATION

Related Product Information for SUMO1 elisa kit	Description: This Mouse Small ubiquitin-related modifier 1 (SUMO1) ELISA Kit employs a two-site sandwich ELISA to quantitate SUMO1 in samples. An antibody specific for SUMO1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any SUMO1 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SUMO1 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SUMO1 bound in the initial step. The color development is stopped and the intensity of the color is measured. Background: Small ubiquitin-related modifier 1 is a member of the SUMO protein family. It functions in a manner similar to ubiquitin in that it is bound to target proteins as part of a post-translational modification system. However, unlike ubiquitin which targets proteins for degradation, this protein is involved in a variety of cellular processes, such as nuclear transport, transcriptional regulation, apoptosis, and protein stability. It is not active until the last four amino acids of the carboxy-terminus have been cleaved off. Several pseudogenes have been reported for this gene. Alternate transcriptional splice variants encoding different isoforms have been characterized. Sequence analysis of the corresponding gene, termed UBL1, revealed that the gene encodes a 101-amino acid polypeptide with homology to ubiquitin and other ubiquitin-like proteins.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	6678489
NCBI GeneID	22218
NCBI Accession #	NP_033486.1
NCBI GenBank Nucleotide #	NP_033486.1
UniProt Accession #	P63166
NCBI Official Full Name	small ubiquitin-related modifier 1
NCBI Official Synonym Full Names	small ubiquitin-like modifier 1
NCBI Official Symbol	Sumo1
NCBI Official Synonym Symbols	GMP1; PIC1; SMT3; Ubl1; SMTP3; Smt3C; SMT3H3; SUMO-1; SENTRIN
NCBI Protein Information	small ubiquitin-related modifier 1
UniProt Protein Name	Small ubiquitin-related modifier 1
UniProt Gene Name	Sumo1
UniProt Synonym Gene Names	Smt3c; Smt3h3; Ubl1; SUMO-1; Smt3C
UniProt Entry Name	SUMO1_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.