

Mouse Antigen Peptide Transporter 2 (TAP2) ELISA Kit | TAP2 elisa kit

Mouse Antigen Peptide Transporter 2 (TAP2) ELISA Kit

CATALOG #
AAA332351

REACTIVITY
Mouse


PRODUCT SYNONYMS

Antigen Peptide Transporter 2 (TAP2); Mouse Antigen Peptide Transporter 2 (TAP2) ELISA Kit; TAP2; DAAP-57C1.2; ABC18; ABCB3; APT2; D6S217E; PSF2; RING11; ABC transporter; MHC 2; ATP-binding cassette; sub-family B (MDR/TAP); member 3; antigen peptide transporter 2; peptide supply factor 2; peptide tr; TAP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	Mouse Antigen peptide transporter 2 (TAP2) ELISA Kit has high sensitivity and excellent specificity for detection of Mouse TAP2. No significant cross-reactivity or interference between Mouse TAP2 and analogues was observed.
Sequence Length	702
Assay Type	Sandwich ELISA (Quantitative)
Samples	Cell Culture Supernatants, Serum, Plasma, Other Biological Fluids
Detection Method	Colorimetric
Usage Notes	* Do not mix components from different kit lots or use reagents beyond the kit expiration date. * Allow all reagents to warm to room temperature for at least 30 minutes before opening. * Pre-rinse the pipet tip with reagent, use fresh pipet tips for each sample, standard and reagent to avoid contamination. * Unused wells must be kept desiccated at 4 degree C in the sealed bag provided. * Mix Thoroughly is very important for the result. It is recommended using low frequency oscillator or slight hand shaking every 10 minutes. * It is recommended that all samples and standards be assayed in duplicate or triplicate.
Preparation and Storage	Store at 2-8 degree C.

IMAGES



Standard Curve (Sample)

Fig.1. Mouse Antigen peptide transporter 2 (TAP2) Standard Curve.

ADDITIONAL INFORMATION

Related Product Information for TAP2 elisa kit	Description: This Mouse Antigen peptide transporter 2 (TAP2) ELISA Kit employs a two-site sandwich ELISA to quantitate TAP2 in samples. An antibody specific for TAP2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TAP2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TAP2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TAP2 bound in the initial step. The color development is stopped and the intensity of the color is measured. Background: The membrane-associated protein encoded by Tap2 is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. The protein encoded by Tap2 is involved in antigen presentation. This protein forms a heterodimer with Tap1 in order to transport peptides from the cytoplasm to the endoplasmic reticulum. Mutations in the human gene may be associated with ankylosing spondylitis, insulin-dependent diabetes mellitus, and celiac disease.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	239915963
NCBI GeneID	21355
NCBI Accession #	NP_035660.3
NCBI GenBank Nucleotide #	NP_035660.3
UniProt Accession #	P36371
NCBI Official Full Name	antigen peptide transporter 2
NCBI Official Synonym Full Names	transporter 2, ATP-binding cassette, sub-family B (MDR/TAP)
NCBI Official Symbol	Tap2
NCBI Official Synonym Symbols	Y1; jas; APT2; Ham2; MTP2; PSF2; ABC18; Abcb3; Ham-2; Tap-2; RING11; AI462429
NCBI Protein Information	antigen peptide transporter 2
UniProt Protein Name	Antigen peptide transporter 2
UniProt Gene Name	Tap2
UniProt Synonym Gene Names	Abcb3; Ham-2; Ham2; APT2
UniProt Entry Name	TAP2_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.