

Micro Tissue Inorganic Phosphorus Assay Kit

Micro Tissue Inorganic Phosphorus Assay Kit

CATALOG #
AAA332625

PRODUCT SYNONYMS

Micro Tissue Inorganic Phosphorus; Micro Tissue Inorganic Phosphorus Assay Kit; Inorganic Phosphorus; Micro Tissue Inorganic Phosphorus assay kit

PRODUCT SPECIFICATIONS

Samples	Animal and Plant Tissue
Usage Note	The Phosphorus Fixative needs to be prepared before use, and is limited to use on the same day, suggest completing the test within 40 min.
Preparation and Storage	Storage at 4 degree C and protected from light upon receipt. Kit has a storage time of 12 months form receipt.

ADDITIONAL INFORMATION

Related Product Information for Micro Tissue Inorganic Phosphorus assay kit	Principle of the Assay: Inorganic phosphorus mainly refers to inorganic phosphate radicals, which participate in various metabolic processes in organisms, including energy metabolism, nucleic acid metabolism, protein phosphorylation and dephosphorylation, and can also promote the synthesis, transformation, and transportation of carbohydrates. Ammonium molybdate reacts with inorganic phosphate to form a substance with a characteristic absorption peak at 660 nm. By measuring the absorbance at 660 nm, the inorganic phosphorus content can be calculated.
Product Categories/Family for Micro Tissue Inorganic Phosphorus assay kit	Cytology; Cell Metabolism; Biochemical Assay Kit

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/assay-kits/micro-tissue-inorganic-phosphorus/332625>