

beta-glucuronidase (beta-GD) Assay Kit | beta-GD assay kit

Micro beta-glucuronidase (beta-GD) Activity Assay Kit

CATALOG #
AAA332666


PRODUCT SYNONYMS

beta-glucuronidase (beta-GD); Micro beta-glucuronidase (beta-GD) Activity Assay Kit; beta-GD assay kit

PRODUCT SPECIFICATIONS

Samples	Animal tissues
Kit Components	Extraction Buffer Reagent I Reagent II Reagent III Reagent IV
Features & Benefits	The verified samples are complete in type and easy to operate.
Usage Notes	If not assayed immediately, samples can be stored at -80 degree C for one month. It is recommended to perform several dilutions of your sample to ensure the readings are within the standard value range. Fresh samples are necessary for good results. If don't perform the assay at the same time, it's better to complete the Sample Preparation step before storing the samples.
Preparation and Storage	Storage at -20 degree C and Keep from light immediately upon receipt. Kit has a storage time of 6 months from receipt. Refer to list of materials supplied for storage conditions of individual components. Shipping: Gel pack with blue ice.

IMAGES



Application Data

ADDITIONAL INFORMATION

Related Product Information for beta-GD assay kit	beta-glucuronidase (beta-GD), a matrix degradation enzyme involved in tumor invasion and metastasis, is widely found in animal tissues. It has physiological functions such as hydrolyzing sterol glucuronic acid and acidic mucopolysaccharide. The content of this enzyme is high in hepatocytes. In addition, it is rich in gastric cancer, so the determination of beta-GD activity in gastric juice is of great significance for the study of gastric cancer. Micro beta-glucuronidase (beta-GD) Activity Assay Kit can detect animal tissues. In this kit, beta-GD catalyzed phenol beta-D-glucuronic acid to produce free phenolphthalein, and the enzyme activity was determined by determining the content of phenol.
Product Categories/Family for beta-GD assay kit	Cell analysis; Cell metabolism detection kit

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/assay-kits/beta-glucuronidase-beta-gd/332666>