

Lipase (LPS) Assay Kit | LPS assay kit

Micro Lipase (LPS) Activity Assay Kit

CATALOG #
AAA332689


PRODUCT SYNONYMS

Lipase (LPS); Micro Lipase (LPS) Activity Assay Kit; LPS assay kit

PRODUCT SPECIFICATIONS

Samples	Animal and Plant Tissues, Cells, Serum (Plasma)
Kit Components	Reagent I Reagent II Reagent III Reagent IV Standard
Features & Benefits	The verified samples are complete in type and easy to operate.
Usage Notes	If not assayed immediately, samples can be stored at -80 degree C for one month. It is recommended to perform several dilutions of your sample to ensure the readings are within the standard value range. Fresh samples are necessary for good results. If don't perform the assay at the same time, it's better to complete the Sample Preparation step before storing the samples.
Preparation and Storage	Storage at 4 degree C and Keep from light immediately upon receipt. Kit has a storage time of 6 months from receipt. Refer to list of materials supplied for storage conditions of individual components. Shipping: Gel pack with blue ice.

IMAGES



Application Data

ADDITIONAL INFORMATION

Related Product Information for LPS assay kit	Lipase (LPS), also known as glycerol ester hydrolase, catalyzes the hydrolysis of triglycerides to produce free fatty acids and glycerol (or di- and mono-glycerides). LPS is widely distributed in various organisms. Elevated levels of LPS in the serum are commonly observed in pancreatitis and pancreatic cancer. Micro Acid Phosphatase (ACP) Activity Assay Kit provides a simple, convenient, and rapid method for detecting LPS activity, applicable to samples such as animal and plant tissues, cells, serum, and plasma. It works by having LPS catalyze the hydrolysis of oil esters into fatty acids, and by measuring the rate of fatty acid production using the copper soap method, the LPS activity can be calculated.
Product Categories/Family for LPS assay kit	Cell analysis ; Cell metabolism detection kit

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/assay-kits/lipase-lps/332689>