

Chymotrypsin Assay Kit

Micro Chymotrypsin Activity Assay Kit

CATALOG #
AAA332697


PRODUCT SYNONYMS

Chymotrypsin; Micro Chymotrypsin Activity Assay Kit; Chymotrypsin assay kit

PRODUCT SPECIFICATIONS

Samples	Animal and plant tissues, bacteria and cells, serum (plasma), or other fluids.
Kit Components	Reagent I Reagent II
Features & Benefits	The verified samples are complete in type and easy to operate.
Usage Notes	If not assayed immediately, samples can be stored at -80 degree C for one month. It is recommended to perform several dilutions of your sample to ensure the readings are within the standard value range. Fresh samples are necessary for good results. If don't perform the assay at the same time, it's better to complete the Sample Preparation step before storing the samples.
Preparation and Storage	Storage at ? degree C and Keep from light immediately upon receipt. Kit has a storage time of 6 months from receipt. Refer to list of materials supplied for storage conditions of individual components. Shipping: Gel pack with blue ice.

IMAGES



Application Data

ADDITIONAL INFORMATION

Related Product Information for Chymotrypsin assay kit	Chymotrypsin, also known as chymotrypsin, is a proteolytic enzyme secreted by the pancreas that can quickly break down denatured proteins. The function of chymotrypsin is similar to that of trypsin, but it has advantages such as strong decompositionability, low toxicity, and minimal adverse reactions. In clinical practice, chymotrypsin is used for sputum dilution and is effective for both purulent and non purulent sputum; It is also used for wound healing after trauma surgery, such as cataract removal. Chymotrypsin catalyzes the hydrolysis of ATEE, and the product exhibits characteristic light absorption at 237 nm; Calculate chymotrypsin activity by measuring the rate of increase in 237 nm light absorption.
Product Categories/Family for Chymotrypsin assay kit	Cell analysis ; Cell metabolism detection kit

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

