

Pepsase Assay Kit

Micro Pepsase Activity Assay Kit

CATALOG #
AAA332698


PRODUCT SYNONYMS

Pepsase; Micro Pepsase Activity Assay Kit; Pepsase assay kit

PRODUCT SPECIFICATIONS

Samples	Animal tissues
Kit Components	Reagent I Reagent II Reagent III Reagent IV Reagent V Reagent VI Standard
Features & Benefits	The verified samples are complete in type and easy to operate.
Usage Notes	If not assayed immediately, samples can be stored at -80 degree C for one month. It is recommended to perform several dilutions of your sample to ensure the readings are within the standard value range. Fresh samples are necessary for good results. If don't perform the assay at the same time, it's better to complete the Sample Preparation step before storing the samples.
Preparation and Storage	Storage at 4 degree C and Keep from light immediately upon receipt. Kit has a storage time of 6 months from receipt. Refer to list of materials supplied for storage conditions of individual components. Shipping: Gel pack with blue ice.

IMAGES



Application Data

ADDITIONAL INFORMATION

Related Product Information for Pepsase assay kit	Pepsin is secreted by the main cells of the gastric mucosa and breaks down proteins in food into small peptides. Generally used for the identification of neuropathic hypoacidosis, chronic gastritis, chronic gastric dilatation, chronic duodenitis and other symptoms will also cause the decrease of pepsin secretion. TM Micro Pepsase Activity Assay Kit can be used to detect biological samples such as animal tissues. In the kit, pepsase can catalyze the hydrolysis of hemoglobin, and the hydrolysate appears blue after reacting with Folin reagent, and the color of pepsin is proportional to the activity of pepsin in a certain range.
Product Categories/Family for Pepsase assay kit	Cell analysis ; Cell metabolism detection kit

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/assay-kits/pepsase/332698>